## English MICRO-MANIPULATION PIPETTES

**NOTE:** One cell mouse embryo tested and passed with 75% or greater blastocyst rate. USP endotoxin (LAL) tested and passed with 20 EU's or less per device. Testing is conducted on a lot-to-lot (batch) basis.

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

**PRECAUTIONS:**

Assisted reproduction microtools (and any other accessories used during these procedures) which contain sperm, zygotes, pre embryos and/or embryos should be comprised of embryo compatible materials.

Chromosomal and congenital abnormalities are complications which may be expected due to the advanced maternal age of the patient population, defects induced through ovulation stimulation and/or manipulation of gametes. When maternal age adjustments are made, the incidence of abnormalities associated with assisted reproduction procedures approximates those of the general population.

**MIKRO-INJECTION PIPETTE (K-MPIP)**

Used for the intracytoplasmic single sperm injection of oocytes. Intended for one-time use.

**CONTRAINDICATIONS:** There are no known contraindications for the use of this device.
**SUGGESTED INSTRUCTION FOR USING MICRO-INJECTION PIPETTE:**

- Preparation For Use:
  - Hold storage container from its ends, between thumb and forefinger.
  - With free hand, open container by lifting lid.
  - Gently hold pipette near its base and lift pipette straight out of holder.
- NOTE: Do not allow pipette to touch the walls or lid of the container during removal. Contact with the container wall or lid could damage the tip of the pipette.**
- Fit micro-injection pipette into tool holder of micromanipulator. Attach to micro-injector.
- Aspirate spermatoza tall-first into tip of micro-injection pipette.
- Once oocyte has been immobilized using the holding pipette, insert tip of micro-injection pipette through the zona pellucida and the oolemma into the ooplasm of the oocyte.
- Inject single sperm along with 1-2 picoliter (pl) of medium into the ooplasm.
- Gently withdraw micro-injection pipette from oocyte.
- Repeat procedure until all oocytes have been injected.
- Remove micro-injection pipette from tool holder and discard.

**-HOLDING PIPETTE (K-HPIP)**

Used to hold an oocyte, embryo, or blastocyst in position with the application of vacuum during intracytoplasmic single sperm injection or assisted hatching/zona drilling. Intended for one-time use.

**CONTRAINDICATIONS:** There are no known contraindications for the use of this device.
**SUGGESTED INSTRUCTION FOR USING HOLDING PIPETTE:**

- Preparation For Use:
  - Hold storage container from its ends, between thumb and forefinger.
  - With free hand, open container by lifting lid.
  - Gently hold pipette near its base and lift pipette straight out of holder.

**NOTE: Do not allow pipette to touch the walls or lid of the container during removal. Contact with the container wall or lid could damage the tip of the pipette.**

- Fit holding pipette into tool holder of micromanipulator. Attach to micro-injector.
- Maneuver Petri dish containing oocytes under microscope until a single oocyte is visible.
- Maneuver tip of holding pipette until it is in position at the selected oocyte. Apply negative pressure to immobilize oocyte.
- Maintain pressure until micro-injection process is complete.
- Release pressure and free oocyte.
- Remove holding pipette from tool holder and discard.

**-ASSISTED HATCHING /ZONA DRILLING PIPETTE (K-AHP)**

Used to make a hole in the zona pellucida to enable embryo assisted hatching. Intended for one-time use.

**CONTRAINDICATIONS:** There are no known contraindications for the use of this device.
**SUGGESTED INSTRUCTION FOR USING ASSISTED HATCHING /ZONA DRILLING PIPETTE:**

- Preparation For Use:
  - Hold storage container from its ends, between thumb and forefinger.
  - With free hand, open container by lifting lid.
  - Gently hold pipette near its base and lift pipette straight out of holder.
- NOTE: Do not allow pipette to touch the walls or lid of the container during removal. Contact with the container wall or lid could damage the tip of the pipette.**
- Remove pipette from tool holder and discard.
- Remove pipette from tool holder and discard.

**-PARTIAL ZONA DISSECTION PIPETTE (K-PZPD)**

Used to create a slit or slits in the zona pellucida to enable assisted hatching and/or cell aspiration.

**PARTIAL ZONA DISSECTION PROCEDURE**

- Pierce zona pellucida with tip of pipette and withdraw pipette.
- Remove pipette from tool holder and discard.

**PRE-IMPLANTATION GENETIC DIAGNOSIS PIPETES**

- Polar Body Biopsy Pipette (K-PBPB)**
- Embryo Biopsy Pipette (K-EBPH)**

Used for the aspiration of cells (i.e. polar bodies, blastomeres, or trophoctoderm cells) for pre-implantation genetic diagnosis. Intended for one-time use.
**NOTE: These tools are indicated for polar body or embryo biopsy, which may be done in order to perform pre-implantation genetic diagnosis (PGD) on the genetic material in the biopsied cells.** Tests for PGD are currently developed, and performance characteristics determined, by individual laboratories for their own use. The performance of these tests may vary depending on the particular assay and disease evaluated. Currently, these tests have not been cleared or approved by the United States Food and Drug Administration.

**SUGGESTED INSTRUCTION FOR USING PRE-IMPLANTATION GENETIC DIAGNOSIS PIPETES:**

- Preparation For Use:
  - Hold storage container from its ends, between thumb and forefinger.
  - With free hand, open container by lifting lid.
  - Gently hold pipette near its base and lift pipette straight out of holder.
- NOTE: Do not allow pipette to touch the walls or lid of the container during removal. Contact with the container wall or lid could damage the tip of the pipette.**
- Fit PGD pipette into tool holder of micromanipulator. Attach to micro-injector.
- Advance the pipette through the existing opening of the zona pellucida and remove the desired cells.

**NOTE: An existing opening is not required when using the Polar Body Biopsy Pipette because it has a beveled spike on the tip to assist in accessing the polar bodies.**


- Remove PGD pipette from tool holder and discard.
- Genetic analysis of cells can now be performed.

- TESTICULAR SPERM EXTRACTION PIPETTE (K-TSEP)**
- Used to remove sperm from a sample of testicular tissue. Intended for one-time use.
**CONTRAINDICATIONS:** There are no known contraindications for the use of this device.
**SUGGESTED INSTRUCTION FOR USING TESTICULAR SPERM EXTRACTION PIPETTE:**

- Obtain testicular tissue sample.
- Tease the tissue away from the sperm with pipette.
- Aspirate sperm and move to new area with pipette.
- Expel sperm ready for ICSI.

**HOW SUPPLIED**

Supplied by gamma irradiation in peel-open pouches. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

North America Representative	
COOK MEDICAL INCORPORATED	Phone: 812 339-2235
P.O. Box 4195	Toll Free: 800 457-4500
Bloomington, IN 47402-4195 USA	Toll Free Fax: 800 554-8335
©COOK 2014	

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