en
20
cs
24
da
28
de
33
el
37
es
42
fi
47
fr
51
hr
56
hu
60
it
65
lt
69
nl
74
no
78
pl
82
pt
87
sl
91
sv
96
tr
100

Frova Intubating Introducer

Instructions for Use

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Frova intubationsindfører

Brugsanvisning

Frova Intubationshilfe

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Kullanma Talimatı

TABLE OF CONTENTS

1. DEVICE DESCRIPTION	20
1.1 Construction of Device.....	20
1.2 Performance Characteristics	20
1.3 Specifications.....	21
1.4 Device Compatibility.....	21
1.5 Patient Population	21
1.6 Intended User	21
1.7 Contact with Body Tissue	21
1.8 Operating Principle.....	21
1.9 Expected Lifetime of the Device.....	21
2. INTENDED USE	21
3. INDICATIONS FOR USE	22
4. CLINICAL BENEFITS	22
5. CONTRAINDICATIONS	22
6. WARNINGS	22
6.1 General Warnings	22
6.2 Sterile and Single Use.....	22
6.3 Malfunctions/Changes in Performance	22
7. PRECAUTIONS	22
7.1 General Precautions	22
7.2 Training Requirements for User	22
7.3 Requirements for Facilities.....	23
8. POTENTIAL ADVERSE EVENTS	23
9. STERILIZATION	23
10. HOW SUPPLIED.....	23
11. INSPECTION OF DEVICE	23
12. INSTRUCTIONS FOR USE.....	23
12.1 Patient Preparation.....	23
12.2 Step-by-Step Guide for Frova Introducer	23
12.3 Step-by-Step Guide for Rapi-Fit Adapters.....	23
13. REQUIRED EQUIPMENT AND ACCESSORIES.....	23
13.1 Materials Required	23
13.2 Materials Recommended	23
14. DISPOSAL OF DEVICE	24
15. REFERENCES	24
16. PATIENT COUNSELING.....	24
17. SERIOUS INCIDENT REPORTING	24

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FROVA INTUBATING INTRODUCER

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

STERILE—DO NOT RESTERILIZE—SINGLE USE ONLY.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

1. DEVICE DESCRIPTION

1.1 Construction of Device

The Frova Intubating Introducer (hereinafter referred to as the Frova Introducer) is a 14 French, 70 cm long radiopaque introducer with centimeter markings and a blunt, curved tip that can be passed into the trachea. The introducer has a through lumen design with two distal sideports to ensure adequate airflow. The Frova Introducer can be supplied with a stiffening stylet and two Rapi-Fit® Adapters for connection to an oxygen source. (Fig. 1) The different configurations are shown in Table 1.

Table 1 – Catalogue Numbers and Components

Catalogue Number	Components
C-CAE-14.0-70-FII	Frova Introducer, a stiffening stylet, a Rapi-Fit Adapter with a 15 mm connector, and a Rapi-Fit Adapter with a female Luer lock connector. Supplied in single units.
C-CAE-14.0-70-FI	Frova Introducer, a Rapi-Fit Adapter with a 15 mm connector, and a Rapi-Fit Adapter with a female Luer lock connector. Supplied in single units.
C-CAE-14.0-70-FIC	Frova Introducer only. Supplied as a box of ten introducers.
C-CAE-14.0-70-FIC-SPOPS	Frova Introducer only. Supplied in single units and packaged in a U-shape for ease of transportation.

1.2 Performance Characteristics

The Frova Introducer is designed to facilitate airway intubation and allows for transient oxygenation and/or jet ventilation during the initial intubation.

The Frova Introducer has two distal sideports, which help ensure sufficient airflow and prevent the Frova Introducer from whipping during jet ventilation. The removable Rapi-Fit Adapters enable the use of an oxygen source if necessary.

The Rapi-Fit Adapters (15 mm or female Luer lock connector) are designed for transient oxygenation and jet ventilation. In Table 2, the ratio of the inspiratory time and expiratory time (I:E ratio), delivered minute volume, and average maximum airway pressure are given for a high-pressure oxygen source tested in adult patients with healthy lung tissue.

Table 2 – High-Pressure Oxygen Source Specifications

Device	I:E Ratio	Average Delivered Minute Volume (L/min)*	Measured Average Maximum Airway Pressure (cm H ₂ O)*
Frova Introducer and Rapi-Fit Adapter with female Luer lock connector	1:4	8.6	12.1**

* The presented numbers are mean values and the following test conditions were used in an active model: Input pressure set to 50 psi, Lung Compliance: 100 mL/cm H₂O, Resistance: 3 cm H₂O/(L/s), tested with ASL 5000, Ingmar Medical, Ltd.
** Measured Average Maximum Airway Pressure did not exceed 14 cm H₂O in the testing performed.

1.3 Specifications

The specifications for the four different Frova configurations are shown in Table 3.

Table 3 – Specifications for Frova Configurations

Catalogue Number	Outer Diameter (Fr/mm)	Length (cm)	Minimum Inner Diameter of Endotracheal Tube (mm)	Stiffening Stylet (cm)	Rapi-Fit w/Female Luer Lock Connector	Rapi-Fit w/15 mm Connector
C-CAE-14.0-70-FII	14.0/4.7	70	6	60	Yes	Yes
C-CAE-14.0-70-FI	14.0/4.7	70	6	No	Yes	Yes
C-CAE-14.0-70-FIC	14.0/4.7	70	6	No	No	No
C-CAE-14.0-70-FIC-SPOPS	14.0/4.7	70	6	No	No	No

1.4 Device Compatibility

The Frova Introducer is compatible with:

- 14 Fr Rapi-Fit Adapter with a 15 mm connector
- 14 Fr Rapi-Fit Adapter with a female Luer lock connector
- Single lumen endotracheal tubes with an inner diameter of 6 mm or larger

Do not use the Frova Introducer with double lumen endotracheal/endobronchial tubes.

1.5 Patient Population

The Frova Introducer, with or without the Rapi-Fit Adapters, can be used in all patients eligible for endotracheal intubation with placement of a single lumen endotracheal tube with an inner diameter of 6 mm or larger and where visualization of the glottis is inadequate e.g., as in grade II and III Cormack-Lehane views.

If the physician evaluates a need for transient oxygenation or jet ventilation during initial intubation, the Rapi-Fit Adapters may be utilized with the Frova Introducer to facilitate this. Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume.

1.6 Intended User

The Frova Introducer has been designed for use by physicians trained and experienced in airway management techniques.

1.7 Contact with Body Tissue

The Frova Introducer is a surface medical device with contact to the mucosal membrane.

1.8 Operating Principle

The Frova Introducer facilitates placement of a single-lumen endotracheal tube. The blunt, curved tip of the Frova Introducer can be passed into the trachea to guide placement of the endotracheal tube, and the centimeter markings facilitate accurate placement of the endotracheal tube.

The stiffening stylet is designed to increase rigidity to the Frova Introducer during placement if preferred by the physician.

If the need for transient oxygenation or jet ventilation occurs, the Rapi-Fit Adapters can be attached to the Frova Introducer to facilitate this in connection with an oxygen source. If an oxygen source is used for insufflation, begin at a lower pressure, and work up gradually. Always make sure that there is enough egression of the insufflated gas volume. Pulse oximetry and oral air flow should be carefully monitored as well. In cases of airway obstruction, gas egression from the patient’s lungs may require more time.

Standard techniques for placement of intubating introducers should be employed.

1.9 Expected Lifetime of the Device

The Frova Introducer is a non-surgical invasive device intended for transient use of less than 60 minutes in the trachea.

2. INTENDED USE

The Frova Intubating Introducer is intended to facilitate endotracheal intubation in patients where the visualization of the glottis is inadequate. The Rapi-Fit Adapters are intended to connect the Frova Intubating Introducer to an oxygen source.

3. INDICATIONS FOR USE

The Frova Intubating Introducer is indicated to facilitate endotracheal intubation in patients where visualization of the glottis is inadequate e.g., as in grade II and III Cormack-Lehane views. This indication applies for introduction of a single lumen endotracheal tube with an inner diameter of 6 mm or larger. The Frova Intubating Introducer is also designed for transient oxygenation or jet ventilation in combination with 14 Fr compatible Rapi-Fit Adapters (15 mm and female Luer lock connector). The Rapi-Fit Adapters should only be used when oxygen requirements are high.

4. CLINICAL BENEFITS

The Frova Introducer facilitates a high overall intubation success rate, safe intubation, and high first-attempt intubation success in patients where visualization of the glottis is inadequate e.g., as in grade II and III Cormack-Lehane views. Furthermore, it can be used with both direct and video laryngoscopy. Transient oxygenation through the Frova Introducer is also possible. When the Frova Introducer is used as intended, the analysis of the clinical data demonstrates that the incidence and type of potential adverse events related to the Frova Introducer are acceptable and that the clinical benefits to the patient population are considered to outweigh the risks.

5. CONTRAINDICATIONS

Do not use the Frova Introducer, if the epiglottis cannot be visualized when performing laryngoscopy, i.e., grade IV Cormack-Lehane laryngoscopy classification.

6. WARNINGS

6.1 General Warnings

- Only use a single lumen endotracheal tube with an inner diameter of 6 mm or larger in combination with the Frova Introducer.
- Do not use the Frova Introducer with double lumen endotracheal/endobronchial tubes.

6.2 Sterile and Single Use

- Do not use the device if the sterile packaging is damaged or unintentionally opened before use. If a contaminated device is introduced to the airways, it may cause an infection and lead to severe harm.
- This single-use device is not designed for reuse. If the device is reused it could pose risks of cross contamination with microbiological agents. Attempts to reprocess (resterilize) and/or to re-use may lead to device failure and/or transmission of disease.

6.3 Malfunctions/Changes in Performance

- Care must be taken when introducing or removing the Frova Introducer from the endotracheal tube as contact between sharp edges on the internal surface of the endotracheal tube may cause small fragments to be shaved off the Frova Introducer during introduction/removal.
- The use of removable Rapi-Fit Adapters permits the utilization of a high and low-pressure oxygen source, if necessary, during the procedure. Use of the Rapi-Fit Adapter for oxygenation or jet ventilation may be associated with a risk of barotrauma.
- Ensure that the Rapi-Fit Adapter is securely connected to the Frova Introducer before proceeding to oxygenation or jet ventilation. Failure to properly secure the adapter to the Frova Introducer may result in reduced oxygen support.
- If a high-pressure oxygen source is used for insufflation (e.g., a jet ventilator), begin at a lower pressure and increase gradually. Rising chest wall, pulse oximetry, and oral air flow should be carefully monitored. **Table 2** provides information on ventilation with a high-pressure oxygen source.

7. PRECAUTIONS

7.1 General Precautions

- Blunt trauma to the thorax may cause injuries distal to the glottis, which may complicate the endotracheal intubation. When possible, imaging of the thorax of trauma patients is recommended prior to intubation to minimize the risk of further injuries to the trachea or main bronchi during intubation.
- Care must be taken not to cause injury to the epiglottis and glottis. Perforation of the sinus pyriformis, trachea, or bronchus must be avoided.
- Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume.
- To avoid barotrauma and/or pneumothorax, examine the patient's anatomy to help determine the optimal placement for the Frova Introducer. Ensure that the Frova Introducer is in the physician-preferred location relative to the carina by referencing its centimeter markings.
- To facilitate oxygenation or jet ventilation through the Frova Introducer, ensure correct connection of the Rapi-Fit Adapter and the oxygen source.

7.2 Training Requirements for User

The physician must be trained and experienced in airway management techniques.

7.3 Requirements for Facilities

The Frova Introducer can be used anywhere intubation may take place e.g., operating room, intensive care unit, emergency department, or pre-hospital.

8. POTENTIAL ADVERSE EVENTS

- Aspiration
- Barotrauma (e.g., pneumothorax or subcutaneous emphysema)
- Bronchial spasm and laryngospasm
- Device fragments in the airway
- Hemodynamic instability
- Hypercapnia (e.g., kink or misplacement of the Frova)
- Hypoxia/hypoxemia (e.g., kink or misplacement of the Frova)
- Respiratory system injury (e.g., pneumothorax, subcutaneous emphysema, hemorrhage, edema, perforation)

9. STERILIZATION

All devices are sterilized using ethylene oxide (EO) gas.

10. HOW SUPPLIED

Supplied in peel-open packages. Keep the device dry and away from sunlight. Do not use after the “use by” (expiration) date printed on the label.

11. INSPECTION OF DEVICE

Inspect the device thoroughly including all levels of the packaging to verify that there is no damage prior to use. Inspect and confirm that the sterile barrier has not been compromised in any way. Inspect and confirm that the device corresponds to the label and the IFU.

Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

12. INSTRUCTIONS FOR USE

12.1 Patient Preparation

The physician should follow institutional protocols and local guidelines. Always ensure adequate pre-oxygenation of the patient if applicable.

12.2 Step-by-step Guide for Frova Introducer

1. Lubricate an appropriately sized single-lumen endotracheal tube (≥ 6 mm) and the distal tip of the Frova Introducer.
2. The endotracheal tube may be preloaded onto the proximal part of the Frova Introducer.
3. Using laryngoscopy, introduce the tip of the Frova Introducer beyond the epiglottis and advance it towards the glottis.
4. Advance the Frova Introducer approximately 2-3 cm into the trachea and remove the stiffening stylet if applicable (only supplied with C-CAE/-FII). If resistance is encountered, do not force the Frova Introducer. Instead, the device should be gently rotated followed by a careful advancing movement.
5. Continue to advance the Frova Introducer no more than 10 cm into the trachea. The optimal distance is based on the size of the patient.
NOTE: The operator should always ensure correct tracheal position of the Frova Introducer and endotracheal tube and never advance the Frova Introducer beyond the carina.
6. Whilst maintaining position of the Frova Introducer, advance the endotracheal tube over the Frova Introducer into the trachea to the appropriate level.
7. Whilst maintaining position of the endotracheal tube, remove the Frova Introducer and the laryngoscope.
8. Secure the endotracheal tube according to the manufacturer’s labeling and confirm its correct position using standard methods e.g., capnography.
NOTE: If necessary, oxygen may be given during the procedure by utilizing the Rapi-Fit Adapters.

12.3 Step-by-step Guide for Rapi-Fit Adapters

1. Attach the unlocked Rapi-Fit Adapter to the Frova Introducer by pushing the white locking collar forward and locking it into position. (**Fig. 2a**)
2. Attach the Rapi-Fit Adapter to an oxygen source.
3. Remove the Rapi-Fit Adapter by pulling the white locking collar back to release and remove from the Frova Introducer. (**Fig. 2b**)

13. REQUIRED EQUIPMENT AND ACCESSORIES

13.1 Materials Required

- Lubricant for airway intubation

13.2 Materials Recommended

- Single-lumen endotracheal tube with an inner diameter of 6 mm or larger
- Laryngoscope
- Oxygenation source
- Pulse oximeter
- Capnograph
- Stethoscope

14. DISPOSAL OF DEVICE

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines. The stiffening stylet contains lead, which should be handled separately.

15. REFERENCES

These instructions for use are based on experience from physicians and/or their published literature. Refer to your local Cook Medical sales representative for information on available literature.

16. PATIENT COUNSELING

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

Furthermore, the patient should be informed about the residual risks when treated with this device and the potential adverse events related to the use of this device.

17. SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, and shall also be reported to the competent authority where the device was used.

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© COOK 2026
William A. Cook Australia Pty Ltd