EN
6

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
1

EN

- Dilator Radiopaque hydrophilically coated 10.0 Fr, 20 cm long
- Femoral preloaded filter introducer with flexible tip
- Tactile Bump
- Coaxial introducer sheath system. Radiopaque introducer sheath, 7.0 Fr, 65 cm long with radiopaque band and radiopaque introducer dilator with 2 radiopaque markers at the distal end
- Check-Flo Valve
- Three-way Stopcock Plastic
- Cook Celest® Platinum Vena Cava Filter (supplied preloaded)
- Hook
- Primary Legs
- Secondary Legs
- Platinum markers
- Anchor

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







10

6

COOK CELECT® PLATINUM VENA CAVA FILTER SET - FOR FEMORAL VEIN APPROACH

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Cook Celect Platinum Filter Set consists of a paramagnetic filter (30 mm diameter, 49 mm long) with platinum markers, preloaded on a femoral filter introducer with a flexible tip, a 7.0 Fr coaxial introducer system (compatible with a .035 inch wire guide) and a hydrophilically coated 10.0 Fr pre-dilator. The introducer dilator has 8 sideports and two radiopaque markers 30 mm apart (end-to-end).

INTENDED USE

The Cook Celect Platinum Filter implant is intended for the prevention of recurrent pulmonary embolism (PE) via placement in the vena cava in the following situations:

- Pulmonary thromboembolism when anticoagulant therapy is contraindicated;
- Failure of anticoagulant therapy in thromboembolic diseases;
- Emergency treatment following massive PE where anticipated benefits of conventional therapy are reduced; and
- Chronic, recurrent PE where anticoagulant therapy has failed or is contraindicated.

The Cook Celect Platinum Filter implant may be retrieved. Please refer to the Instructions for Use provided with the Günther Tulip Vena Cava Retrieval Set (not included in the filter set).

The product is intended for percutaneous placement via a femoral vein for filtration of inferior vena cava (IVC) blood to prevent PE.

CONTRAINDICATIONS

Filter Placement

- Megacava (diameter of the IVC > 30 mm).
- Diameter of the IVC < 15 mm.
- Extensive thrombus in the vein chosen for approach.
- Vena Cava filters should not be implanted in patients with risk of septic embolism due to the risk of infection. The decision should be based on the patient's individual risk/benefit profile.

Optional Filter Retrieval

- Retrieval of the filter with significant amounts of trapped thrombus (greater than 25% of the volume of the cone).
- Retrieval of the filter for patients with an on-going high risk for PE.

WARNINGS

Filter Placement

- Manipulation of products requires imaging control.
- When injecting contrast media, do not exceed maximum pressure rating of 1000 psi and flowrate of 20 ml/sec. Hand injection is also possible.
- **When placing the premounted filter via the femoral approach into the Check-Flo® valve of introducer sheath, hold the introducer with flexible tip near the end, close to the filter, to avoid kinking of the flexible tip.**
- Do not rotate the preloaded filter inside the introducer system.
- **Excessive force should not be exerted to place filter.** If severe resistance is met when advancing the wire guide, then retract the wire guide and choose a different approach.

Optional Filter Retrieval

- **Excessive force should not be exerted to retrieve the filter.**
- An inferior vena caval evaluation for residual captured thrombus should be performed prior to attempted retrieval.
- Available data from retrievals in a prospective, multicenter study demonstrate that the device can be safely retrieved. Please refer to the **Clinical Studies** section of this booklet for clinical study references to the retrieval of this filter.

PRECAUTIONS

- Possible allergic reactions should be considered.
- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

Filter Placement

- For placement of the filter, the right femoral vein is usually preferred due to the route to the vena cava. The left femoral vein can be used, but is more tortuous. Prior to choosing an approach, assess

the patient's size, anatomy and the location of the venous thrombosis.

- The filter for femoral vein approach is supplied preloaded on the filter introducer. Do not separate the preloaded filter introducer assembly to view or examine the components.
- Do not attempt to reload the filter onto the filter introducer. Any attempt to do so may damage the introducer or the filter.
- Before any injection of contrast media through the sheath, ensure that the sheath hub and handle are connected, or confirm by diagnostic imaging of the sheath tip (or radiopaque marker).
- Once the metal mount point is past the tip of the sheath, the secondary legs of the filter are expanded. The filter may be repositioned only by advancing the filter; retracting the filter could damage the secondary legs or caval wall.

Optional Filter Retrieval

- For filter retrieval, a right jugular vein approach is preferable. An approach via the left jugular vein is possible; however no data are available demonstrating the safety or effectiveness of filter retrieval via the left jugular vein.
- The filter has been designed to be retrieved with the Günther Tulip Vena Cava Filter Retrieval Set, (not included in the filter set). Cook has not performed testing to evaluate the safety or effectiveness of filter retrieval using other retrieval systems.
- Never re-deploy a retrieved filter.
- Filter retrieval is a patient-specific, clinically complex decision; the decision to remove a filter should be based on the patient's individual risk/benefit profile. Retrieve the filter when feasible and clinically indicated.
- The clinical study data suggest that the optimal retrieval period is between 0 and 52 weeks, during which time the probability of successfully retrieving the filter is greater than 90.0%. Filter retrieval is a patient-specific, clinically complex decision; a filter should only be retrieved when feasible and clinically indicated. For all IVC filter types, retrieval becomes more challenging with time. This is due to encapsulation of the filter legs or hook (in a tilted filter) by tissue in-growth. Please refer to the "Clinical Studies" section of this booklet for further details related to filter retrieval and the reference to the clinical study describing the retrievability of this filter.

MRI INFORMATION

Non-clinical testing has demonstrated that the Cook Celect Filter is **MR Conditional** according to ASTM F2503. The MR safety testing conducted on the Cook Celect filter is considered sufficient to support the MR safety of the Celect Platinum filter. A patient with this filter can be scanned safely after placement under the following conditions.

- Static magnetic field of 3.0 Tesla or 1.5 Tesla.
- Maximum spatial magnetic gradient of 1600 Gauss/cm (16.0 T/m) or less.
- Normal operating mode.
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning or less (i.e., per scanning sequence).

Static Magnetic Field

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 Tesla Temperature Rise

In non-clinical testing, the Cook Celect Filter produced a maximum temperature change of 3.7 °C during 15 minutes of MRI (i.e., for one scanning sequence) performed in a 1.5 Tesla MR System (1.5 Tesla/64 MHz, Siemens Magnetom Avanto, NUMARIS/4 syngo MR B17 DHHS) at an MR system reported whole body averaged SAR of 2.05 W/kg (associated with a calorimetry measured whole body averaged value of 1.75 W/kg).

3.0 Tesla Temperature Rise

In non-clinical testing, the Cook Celect Filter produced a maximum temperature change of 4.2 °C during 15 minutes of MRI (i.e., for one scanning sequence) performed in a 3.0 Tesla MR System (3.0 Tesla/128 MHz, Siemens Magnetom Trio, A Tim System, NUMARIS/4 syngo MR B17 DHHS) at an MR system reported whole body averaged SAR of 2.08 W/kg (associated with a calorimetry measured whole body averaged value of 1.71 W/kg).

Image Artifacts

MR image quality may be compromised if the area of interest is within approximately 21 mm of the position of the Cook Celect Filter as found during non-clinical testing using T1-weighted, spin echo and gradient

echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic filter.

POTENTIAL ADVERSE EVENTS

- Damage to the vena cava
- Pulmonary embolism
- Filter embolization
- Vena cava perforation/penetration
- Vena cava occlusion or thrombosis
- Hemorrhage
- Hematoma at vascular access site
- Infection at vascular access site
- Cardiac tamponade
- Filter malpositioning
- Postphlebotic syndrome
- Death

CLINICAL STUDIES


Previously published clinical studies for the Celect filter suggest probable clinical results for the successful retrieval of the Celect Platinum filter.

Sixty-one (61) men and 34 women (mean age: 51 ± 19 years) were enrolled in a prospective, multi-center

registry study to assess the safety of retrieving the Cook Celect Vena Cava Filter in patients with a high risk of pulmonary thromboembolism (PE).

All 95 filters were successfully placed. Two patient deaths during follow-up (one within 24 hours of the procedure and one associated with a recurrent PE) were adjudicated as potentially device-related by an independent Clinical Events Committee. Filter retrieval was attempted when deemed clinically appropriate by the physician. Retrieval was attempted in 58 patients. Imaging at retrieval showed no evidence of IVC perforation and 21 cases of filter strut incorporation into the IVC wall. The mean time to retrieval was 179 days (range 5 - 466 days). Fifty-six (56) retrievals (96.6%) were successful. The two unsuccessful retrievals were attributed to an inability to capture the filter hook due either to excessive filter tilt or to tissue growth causing the hook to become embedded in the endothelium.

A Kaplan-Meier product limit estimate (see Figures below) indicates that the estimated probability of successful retrieval of the Celect filter remains at 100% at up to 50 weeks post-implant and above 74% at 55 weeks post-implant.

Filter Indwell Time Weeks	Kaplan-Meier Estimated Probability of Successfully Retrieving the Celect Filter	Standard Error
4	100%	0.00
12	100%	0.00
26	100%	0.00
52	90%	0.09
60	75%	0.16

Study Publication:

Lyon SM, Riojas GE, Ueberoi R, Patel J, Lipp ME, Plant GR, De Gregorio MA, Günther RW, Voorhees WD, and McCann-Brown JA. Short and long-term retrievability of the Celect Vena Cava Filter: results from a multi-institutional registry. J Vasc Interv Radiol 2009; 20:1441-1448.

INSTRUCTIONS FOR USE

Femoral Approach

General

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed. It is assumed that the operator will use local anesthesia, sedation and analgesia as required.

Preparation

1. Flush the introducer sheath and dilator.
2. Advance the introducer dilator through the middle of the Check-Flo valve on the introducer sheath. Secure the introducer dilator to the introducer sheath by twisting the dilator hub clockwise (**Fig. 2**).
3. Remove the filter protection tube (**Fig. 3**).

Filter Placement

4. Access the chosen femoral vein using the Seldinger technique.
5. If the IVC anatomy is known from recent imaging, then proceed directly to filter placement.
6. If the IVC anatomy (i.e., dimensions and morphology) is not known, perform imaging to confirm a single IVC, measure the IVC diameter, check for thrombus, and establish the position of the renal veins.
7. Place a supportive .035 inch wire guide in the IVC.
8. If necessary, dilate the puncture site with the 10.0 Fr pre-dilator.
9. Remove the dilator and advance the coaxial introducer sheath system over the wire guide until the tip of the sheath lies just caudal to the renal veins.

10. Remove the wire guide.

11. Using either power or hand injection, perform diagnostic imaging to verify position below (caudal to) the renal veins.

WARNING: When using a power injector, do not exceed maximum pressure rating of 1000 psi and flow rate of 20 ml/sec.

WARNING: Before injecting contrast media through the introducer dilator, ensure that the sheath hub and handle are correctly connected, or confirm by diagnostic imaging of the sheath tip (or radiopaque marker).

12. When correct position is established, remove the introducer dilator by twisting the dilator hub counterclockwise (**Fig. 4**).
13. When placing the premounted filter via the femoral approach into the Check-Flo valve of introducer sheath, hold delivery catheter near the end, close to the filter, to avoid kinking of the flexible tip (**Fig. 5**).

WARNING: Do not rotate the preloaded filter inside the introducer system.

14. Advance the filter introducer until the Check-Flo valve contacts the tactile bump. This will place the hook of the filter inside the sheath at the radiopaque band. Verify that the position of the hook is inside the introducer sheath and still below the renal veins (**Fig. 6**).

WARNING: Excessive force should not be exerted to advance the filter through the delivery system.

15. Stabilize the introducer, and withdraw the introducer sheath and connect it to the handle. At this point the filter is expanded, still connected to the filter introducer (**Figs. 7 and 8**).

WARNING: Repositioning is possible only by advancing the filter; retracting the filter could damage the secondary legs or the caval wall.

16. Proper position can now be verified by diagnostic imaging.

WARNING: Do not rotate the released filter inside the vena cava. Doing so may compromise the performance of the filter.

CAUTION: Filter release or injection of contrast medium must not be made unless the metal mounting point is completely free of the sheath. Use the radiopaque band for positioning.

17. Verify that the sheath hub and handle are connected to ensure that the metal mounting point is completely free of the sheath before filter release.
18. When the filter position is correct, push the red safety button to prepare filter release (**Fig. 9**).
19. To release the filter, push the release button completely to ensure proper release of filter. Repositioning of the filter is no longer possible. The filter is now released (**Fig. 10**).
20. Perform diagnostic imaging to verify filter position, and then withdraw the introducer sheath.

NOTE: Hospital standard of care should be followed for removing the introducer sheath and providing hemostasis to prevent bleeding at the vascular access site.

Optional Retrieval Procedure

NOTE: The Cook Select Platinum Filter implant may be retrieved. Please refer to the Instructions for Use provided with the Günther Tulip Vena Cava Retrieval Set (not included in the filter set).

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2014