en 29	Zenith® Fenestrated AAA Endovascular Graft with the H&L-B One-Shot™ Introduction System Instructions for use
cs 38	Endovaskulární graft Zenith® Fenestrated AAA se zaváděcím systémem H&L-B One-Shot™ Návod k použití
da 47	Zenith® Fenestrated AAA endovaskulær protese med H&L-B One-Shot™ indføringssystem Brugsanvisning
de 57	Zenith® Fenestrated endovaskuläre AAA-Prothese mit dem H&L-B One-Shot™ Einführsystem Gebrauchsanweisung
el 68	Ενδοαγγειακό μόσχευμα AAA Zenith® Fenestrated με το σύστημα εισαγωγής H&L-B One-Shot™ Οδηγίες χρήσης
es 80	Endoprótesis vascular para AAA Zenith® Fenestrated con el sistema de introducción H&L-B One-Shot™ Instrucciones de uso
et 90	AAA endovaskulaarne protees Zenith® Fenestrated sisestussüsteemiga H&L-B One-Shot™ Kasutusjuhend
fi 100	Zenith® Fenestrated AAA -endovaskulaarinen siirre ja H&L-B One-Shot™-sisäänvientijärjestelmä Käyttöohjeet
fr 110	Endoprothèse vasculaire Zenith® Fenestrated AAA à système d'introduction H&L-B One-Shot™ Mode d'emploi
hr 120	Endovaskularni transplantat AAA Zenith® Fenestrated sa sustavom za uvođenje H&L-B One-Shot™ Upute za uporabu
hu 130	Zenith® Fenestrated AAA endovaszkuláris graft H&L-B One-Shot™ felvezetőrendszerrel Használati utasítás
it 140	Endoprotesi addominale Zenith® Fenestrated AAA con sistema di introduzione H&L-B One-Shot™ Istruzioni per l'uso
lt 150	„ Zenith® Fenestrated “ AAA endovaskulinis protezas su „H&L-B One-Shot™“ įvedimo sistema Naudojimo instrukcija
lv 160	Zenith® Fenestrated AAA endovaskulārā protēze ar H&L-B One-Shot™ ievadišanas sistēmu Lietošanas pamācība
nl 170	Zenith® Fenestrated AAA endovasculaire prothese met het H&L-B One-Shot™ introductiesysteem Gebruiksaanwijzing
no 180	Zenith® Fenestrated AAA endovaskulært implantat med H&L-B One-Shot™ innføringssystem Bruksanvisning
pl 189	Stent-graft wewnątrznaczyniowy Zenith® Fenestrated AAA z systemem wprowadzającym H&L-B One-Shot™ Instrukcja użycia
pt 200	Prótese endovascular Zenith® Fenestrated AAA com o sistema de introdução H&L-B One-Shot™ Instruções de utilização
sk 210	Endovaskulárny štep Zenith® Fenestrated AAA so zavádzacím systémom H&L-B One-Shot™ Návod na použitie
sl 220	Endovaskularni vsadek Zenith® Fenestrated AAA z uvajalnim sistemom H&L-B One-Shot™ Navodila za uporabo
sv 230	Zenith® Fenestrated AAA endovaskulärt graft med H&L-B One-Shot™ införingssystem Bruksanvisning

TABLE OF CONTENTS

Illustrations	13–28
1. DEVICE DESCRIPTION	29
1.1 Construction of device	29
1.2 Performance characteristics	29
1.2.1 Graft	29
1.2.2 Proximal body graft	29
1.2.3 Distal bifurcated body graft	29
1.2.4 Delivery system	29
1.3 Device compatibility	29
1.3.1 Device sizing options	29
1.3.2 Device diameter sizing guidelines	30
1.3.3 Iliac leg	30
1.3.4 Bridging stents	31
1.3.5 Materials required	31
1.3.6 Materials recommended	31
1.4 Qualitative and quantitative information	31
1.5 Patient population	31
1.6 Intended user	31
1.7 Contact with body tissue	31
1.8 Operating principle	31
1.9 Expected lifetime of device	31
2. INTENDED USE	31
3. INDICATIONS FOR USE	31
4. CLINICAL BENEFITS	31
5. CONTRAINDICATIONS	31
6. WARNINGS	31
6.1 General	31
6.2 Sterile and single use	32
6.3 Patient selection, treatment and follow-up	32
6.4 Implant procedure	32
6.5 Molding balloon use	32
7. PRECAUTIONS	32
7.1 Training requirements for user	32
8. MRI SAFETY INFORMATION	33
9. POTENTIAL ADVERSE EVENTS	33
10. HOW SUPPLIED	33
11. INSPECTION OF THE DEVICE	33
12. DEVICE PREPARATION	33
13. INSTRUCTIONS FOR USE	33
13.1 General use information	33
13.2 Pre-implant review	33
13.3 Patient preparation	33
13.4 Fenestrated system	33
13.4.1 Proximal body preparation/flush	33
13.4.2 Distal bifurcated body graft preparation/flush	33
13.4.3 Vascular access and angiography	34
13.4.4 Proximal body placement	34
13.4.5 Docking of top cap	34
13.4.6 Bridging stent placement and deployment	34
13.4.7 Distal bifurcated body placement	34
13.4.8 Contralateral iliac wire guide placement	35
13.4.9 Distal bifurcated body deployment	35
13.4.10 Iliac leg (contralateral) placement	35
13.4.11 Distal bifurcated body deployment	35
13.4.12 Molding balloon insertion	35
13.5 Final angiogram	35
13.6 Trigger-wire release troubleshooting	35
13.6.1 Alternate proximal body deployment	35
13.6.2 Docking of top cap	36
13.7 Suprarenal stent deployment troubleshooting	36
13.7.1 Proximal body placement with distal attachment	36
13.7.2 Proximal body placement without distal attachment	36
14. IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP	36
14.1 General	36
14.2 Contrast and non-contrast CT recommendations	37
14.3 Abdominal radiographs	37
14.4 Ultrasound	37
14.5 Additional surveillance and treatment	37
15. DISPOSAL OF DEVICES	37
16. PATIENT COUNSELLING INFORMATION	37
17. SERIOUS INCIDENT REPORTING	37

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
Figure 68


Figure 69

ZENITH® FENESTRATED AAA ENDOVASCULAR GRAFT WITH THE H&L-B ONE-SHOT™ INTRODUCTION SYSTEM

Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

1. DEVICE DESCRIPTION

1.1 Construction of device

The Zenith® Fenestrated AAA Endovascular Graft is a modular system consisting of both a fenestrated proximal body and a bifurcated distal body and used in combination with an iliac leg graft. **(Figure 1)**

The grafts are constructed of full-thickness woven polyester fabric sewn to self-expanding stainless steel Cook-Z stents™ with braided polyester and monofilament polypropylene suture.

The Zenith Fenestrated AAA Endovascular Graft is shipped preloaded onto the H&L-B One-Shot Introduction System. Each individual device has its own separate delivery system. The Proximal Component uses a 6.7 mm I.D. (20 French) or 7.3 mm I.D. (22 French) H&L-B One-Shot Introduction System and the Distal component uses a 6.7 mm I.D. (20 French) H&L-B One-Shot Introduction System.

1.2 Performance characteristics

1.2.1 Graft

The modular nature of the Zenith Fenestrated AAA Endovascular Graft (Proximal Body Graft and Distal Bifurcated body graft used in combination with an iliac leg) allows each graft to be more accurately sized to the patient's anatomy.

The grafts are fully stented to provide integrity, stability and the radial force necessary to open the lumen of the graft during deployment and seal the graft to the vessel wall. The combination of the self-expanding Stainless Steel Cook Z-Stents sewn onto the polyester fabric results in exclusion of blood flow to the aneurysm and through the graft lumen, without causing damage to the vessel wall.

1.2.2 Proximal body graft

The bare suprarenal stent at the proximal end assists the graft to remain in position in the vasculature using stainless steel barbs placed at 3 mm circumferential intervals.

Up to three precisely located holes (fenestration(s)), and cut-outs from the proximal margin (scallop(s)) of the graft material allow for patient specific precision treatment as well as allow the proximal margin of the device to sit higher than standard AAA devices and allow uninterrupted blood flow to branch vessels of the aorta such as the renal and superior mesenteric arteries. **(Figure 2)**

Radiopaque markers are positioned at key positions on the graft to facilitate advancement, deployment and withdrawal under fluoroscopy:

- Four in a circumferential orientation within 1 mm of the most superior aspect of the graft material.
- One on the lateral aspect of the most distal stent.
- On the anterior and posterior surfaces of the graft below the level of the fenestrations (The vertical anterior markers, and the horizontal posterior markers should form a cross when correctly oriented).
- Anteriorly in the shape of a tick or check marker.
- Along the edges of the scallops and fenestrations.

1.2.3 Distal bifurcated body graft

Radiopaque markers are positioned at key positions on the graft to facilitate advancement, deployment and withdrawal under fluoroscopy:

- One radiopaque marker at the graft bifurcation.
- One contralateral limb radiopaque marker at the distal end of the contralateral limb. **(Figure 3)**

1.2.4 Delivery system

The Zenith Fenestrated Proximal **(Figure 4)** and Distal **(Figure 5)**

Component are preloaded onto separate H&L-B One-Shot Introduction systems. The delivery system has a sequential deployment method with built-in features to provide continuous control of the graft throughout the deployment procedure.

The Proximal segment (bare suprarenal stent on the Proximal component) and distal segment of the grafts are attached to the delivery system and held by internal wires enabling precise positioning and allowing readjustment of the graft position before full deployment of the graft.

The bare suprarenal stent of the proximal component is constrained within a top cap to further support the precise positioning and readjustment of the graft.

The Zenith Fenestrated Proximal Graft is reduced in diameter by an independent wire tied to diameter reducing ties to allow the graft to be manipulated within the aorta to allow accurate positioning of the graft and enabling the fenestration(s) to line up with the desired arteries.

For added hemostasis, the Captor® Hemostatic Valve can be opened or closed for the introduction and/or removal of ancillary devices into and out of the sheath.

The proximal and distal body graft delivery system features a Flexor® introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.3 Device compatibility

1.3.1 Device sizing options

The Zenith Fenestrated AAA Endovascular Graft and delivery systems are available in the lengths and diameters detailed in tables 1.1 to 1.2.

Reorder number	1.1 Proximal body graft				
	Proximal body diameter	Proximal body length	Introduction sheath		Introduction sheath length
	mm	mm	French size	(I.D. / O.D.) mm	cm
ZFEN-P-1-24-76	24	76	20	6.7 / 7.7	50
ZFEN-P-1-26-76	26	76	20	6.7 / 7.7	50
ZFEN-P-1-28-76	28	76	20	6.7 / 7.7	50
ZFEN-P-1-30-76	30	76	20	6.7 / 7.7	50
ZFEN-P-1-32-76	32	76	20	6.7 / 7.7	50
ZFEN-P-1-24-91	24	91	20	6.7 / 7.7	50
ZFEN-P-1-26-91	26	91	20	6.7 / 7.7	50
ZFEN-P-1-28-91	28	91	20	6.7 / 7.7	50
ZFEN-P-1-30-91	30	91	20	6.7 / 7.7	50
ZFEN-P-1-32-91	32	91	20	6.7 / 7.7	50
ZFEN-P-1-24-106	24	106	20	6.7 / 7.7	50
ZFEN-P-1-26-106	26	106	20	6.7 / 7.7	50
ZFEN-P-1-28-106	28	106	20	6.7 / 7.7	50
ZFEN-P-1-30-106	30	106	20	6.7 / 7.7	50
ZFEN-P-1-32-106	32	106	20	6.7 / 7.7	50
ZFEN-P-1-24-121	24	121	20	6.7 / 7.7	50
ZFEN-P-1-26-121	26	121	20	6.7 / 7.7	50
ZFEN-P-1-28-121	28	121	20	6.7 / 7.7	50
ZFEN-P-1-30-121	30	121	20	6.7 / 7.7	50
ZFEN-P-1-32-121	32	121	20	6.7 / 7.7	50
ZFEN-P-1-34-84	34	84	20	6.7 / 7.7	50
ZFEN-P-1-36-84	36	84	22	7.3 / 8.5	50
ZFEN-P-1-34-99	34	99	20	6.7 / 7.7	50
ZFEN-P-1-36-99	36	99	22	7.3 / 8.5	50
ZFEN-P-1-34-114	34	114	20	6.7 / 7.7	50
ZFEN-P-1-36-114	36	114	22	7.3 / 8.5	50
ZFEN-P-1-34-129	34	129	20	6.7 / 7.7	50
ZFEN-P-1-36-129	36	129	22	7.3 / 8.5	50
ZFEN-P-2-24-94	24	94	20	6.7 / 7.7	50
ZFEN-P-2-26-94	26	94	20	6.7 / 7.7	50
ZFEN-P-2-28-94	28	94	20	6.7 / 7.7	50
ZFEN-P-2-30-94	30	94	20	6.7 / 7.7	50
ZFEN-P-2-32-94	32	94	20	6.7 / 7.7	50
ZFEN-P-2-24-109	24	109	20	6.7 / 7.7	50
ZFEN-P-2-26-109	26	109	20	6.7 / 7.7	50
ZFEN-P-2-28-109	28	109	20	6.7 / 7.7	50
ZFEN-P-2-30-109	30	109	20	6.7 / 7.7	50
ZFEN-P-2-32-109	32	109	20	6.7 / 7.7	50
ZFEN-P-2-24-124	24	124	20	6.7 / 7.7	50
ZFEN-P-2-26-124	26	124	20	6.7 / 7.7	50
ZFEN-P-2-28-124	28	124	20	6.7 / 7.7	50
ZFEN-P-2-30-124	30	124	20	6.7 / 7.7	50
ZFEN-P-2-32-124	32	124	20	6.7 / 7.7	50
ZFEN-P-2-34-107	34	107	20	6.7 / 7.7	50
ZFEN-P-2-36-107	36	107	22	7.3 / 8.5	50
ZFEN-P-2-34-122	34	122	20	6.7 / 7.7	50
ZFEN-P-2-36-122	36	122	22	7.3 / 8.5	50
ZFEN-P-2-34-137	34	137	20	6.7 / 7.7	50
ZFEN-P-2-36-137	36	137	22	7.3 / 8.5	50

Reorder number	Ipsilateral leg diameter mm	Distal body length mm	1.2 Distal bifurcated body graft			
			Distal leg length mm	French size Fr	Introduction sheath (I.D. / O.D.) mm	Introduction sheath length cm
ZFEN-D-12-28-76	12	76	28	20	6.7 / 7.7	40
ZFEN-D-16-28-76	16	76	28	20	6.7 / 7.7	40
ZFEN-D-20-28-76	20	76	28	20	6.7 / 7.7	40
ZFEN-D-12-28-94	12	94	28	20	6.7 / 7.7	40
ZFEN-D-16-28-94	16	94	28	20	6.7 / 7.7	40
ZFEN-D-20-28-94	20	94	28	20	6.7 / 7.7	40
ZFEN-D-12-28-109	12	109	28	20	6.7 / 7.7	40
ZFEN-D-16-28-109	16	109	28	20	6.7 / 7.7	40
ZFEN-D-20-28-109	20	109	28	20	6.7 / 7.7	40
ZFEN-D-12-28-124	12	124	28	20	6.7 / 7.7	40
ZFEN-D-16-28-124	16	124	28	20	6.7 / 7.7	40
ZFEN-D-20-28-124	20	124	28	20	6.7 / 7.7	40
ZFEN-D-12-45-76	12	76	45	20	6.7 / 7.7	40
ZFEN-D-16-45-76	16	76	45	20	6.7 / 7.7	40
ZFEN-D-20-45-76	20	76	45	20	6.7 / 7.7	40
ZFEN-D-24-45-76	24	76	45	20	6.7 / 7.7	40
ZFEN-D-12-45-94	12	94	45	20	6.7 / 7.7	40
ZFEN-D-16-45-94	16	94	45	20	6.7 / 7.7	40
ZFEN-D-20-45-94	20	94	45	20	6.7 / 7.7	40
ZFEN-D-24-45-94	24	94	45	20	6.7 / 7.7	40
ZFEN-D-12-45-109	12	109	45	20	6.7 / 7.7	40
ZFEN-D-16-45-109	16	109	45	20	6.7 / 7.7	40
ZFEN-D-20-45-109	20	109	45	20	6.7 / 7.7	40
ZFEN-D-24-45-109	24	109	45	20	6.7 / 7.7	40
ZFEN-D-12-45-124	12	124	45	20	6.7 / 7.7	40
ZFEN-D-16-45-124	16	124	45	20	6.7 / 7.7	40
ZFEN-D-20-45-124	20	124	45	20	6.7 / 7.7	40
ZFEN-D-24-45-124	24	124	45	20	6.7 / 7.7	40
ZFEN-D-12-62-76	12	76	62	20	6.7 / 7.7	40
ZFEN-D-16-62-76	16	76	62	20	6.7 / 7.7	40
ZFEN-D-20-62-76	20	76	62	20	6.7 / 7.7	40
ZFEN-D-24-62-76	24	76	62	20	6.7 / 7.7	40
ZFEN-D-12-62-94	12	94	62	20	6.7 / 7.7	40
ZFEN-D-16-62-94	16	94	62	20	6.7 / 7.7	40
ZFEN-D-20-62-94	20	94	62	20	6.7 / 7.7	40
ZFEN-D-24-62-94	24	94	62	20	6.7 / 7.7	40
ZFEN-D-12-62-109	12	109	62	20	6.7 / 7.7	40
ZFEN-D-16-62-109	16	109	62	20	6.7 / 7.7	40
ZFEN-D-20-62-109	20	109	62	20	6.7 / 7.7	40
ZFEN-D-24-62-109	24	109	62	20	6.7 / 7.7	40
ZFEN-D-12-62-124	12	124	62	20	6.7 / 7.7	40
ZFEN-D-16-62-124	16	124	62	20	6.7 / 7.7	40
ZFEN-D-20-62-124	20	124	62	20	6.7 / 7.7	40
ZFEN-D-24-62-124	24	124	62	20	6.7 / 7.7	40

1.3.2 Device diameter sizing guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow.

Intended aortic vessel diameter mm	Main body diameter mm	1.3 Proximal body graft diameter sizing	
		French size Fr	Introduction sheath (I.D. / O.D.) mm
19	24	20	6.7 / 7.7
20	24	20	6.7 / 7.7
21	24	20	6.7 / 7.7
21	26	20	6.7 / 7.7
22	26	20	6.7 / 7.7
23	28	20	6.7 / 7.7
24	28	20	6.7 / 7.7
24	30	20	6.7 / 7.7
25	30	20	6.7 / 7.7
26	30	20	6.7 / 7.7
26	32	20	6.7 / 7.7
27	32	20	6.7 / 7.7
28	32	20	6.7 / 7.7
28	34	20	6.7 / 7.7
29	34	20	6.7 / 7.7
29	36	22	7.3 / 8.5
30	36	22	7.3 / 8.5
31	36	22	7.3 / 8.5

Intended iliac vessel diameter mm	Ipsilateral fenestrated leg mm	1.4 Distal ipsilateral graft diameter sizing	
		French size Fr	Introduction sheath (I.D. / O.D.) mm
9	12	20	6.7 / 7.7
10	12	20	6.7 / 7.7
11	12	20	6.7 / 7.7
12	16	20	6.7 / 7.7
13	16	20	6.7 / 7.7
14	16	20	6.7 / 7.7
15	16	20	6.7 / 7.7
15	20	20	6.7 / 7.7
16	20	20	6.7 / 7.7
17	20	20	6.7 / 7.7
18	20	20	6.7 / 7.7
18	24	20	6.7 / 7.7
19	24	20	6.7 / 7.7
20	24	20	6.7 / 7.7
21	24	20	6.7 / 7.7

1.3.3 Iliac leg

The Zenith Fenestrated AAA Endovascular Graft is used in combination with the Zenith® Spiral-Z® AAA Iliac Leg (ZSLE) (Figure 6). Refer to instructions for use for this device.

1.3.4 Bridging stents

It is recommended that in the event that small fenestrations are being utilized, bridging stents may be placed to secure positive alignment. In cases where the physician elects to place a bridging stent within the small fenestration, a balloon-expandable covered stent must be used. Balloon-expandable covered stents, used as bridging stents, shall be indicated for placement within renal/visceral arteries. The length of the bridging stent to be used is dependent on the patient's anatomy. Consult the balloon-expandable covered stents manufacturer's instructions for selecting the appropriately sized balloon-expandable covered stent.

1.3.5 Materials required

The device is required to be used in conjunction with the following. For information on the use of these products, refer to the individual product's instructions for use.

(Not included in modular system):

- Fluoroscope with digital angiography capabilities
- Contrast media
- Syringe
- Heparinised saline solution
- Sterile gauze pads

1.3.6 Materials recommended

The device is recommended to be used in conjunction with the following. For information on the use of these products, refer to the individual product's instructions for use.

(Not included in modular system):

- 0.035 inch (0.89 mm) extra stiff wire guide, 260 cm; for example:
 - Cook Medical's Amplatz Ultra Stiff Wire Guide (AUS)
 - Cook Medical's Lunderquist® Extra-Stiff Wire Guide (LES)
- 0.035 inch (0.89 mm) standard soft-tipped wire guide; for example:
 - Cook Medical's 0.035 inch soft-tipped Rosen wire guides
- Molding Balloons; for example:
 - Cook Medical's Coda® Balloon Catheters
- Introducer sets; for example:
 - Cook Medical's Check-Flo® Introducer Sets
 - Cook Medical's Flexor Guiding Sheaths
- Sizing catheter; for example:
 - Cook Medical's Aorous™ Centimeter sizing catheter
- Angiographic radiopaque tip catheters; for example:
 - Cook Medical's angiographic catheters
 - Cook Medical's Royal Flush® catheters
- 18 gauge or 19 ultra-thin gauge entry needles; for example:
 - Cook Medical's one part percutaneous needles
- Balloon-expandable covered bridging stents

1.4 Qualitative and quantitative information

		Weight (grams)	
Implant material		ZFEN-P	ZFEN-D
Stainless steel	Stent	Up to 4.26	Up to 3.21
	Graft fabric	Up to 2.63	Up to 2.53
Polyester	Suture	Up to 0.23	Up to 0.22
	Suture	Up to 0.28	Up to 0.22
Polypropylene	Suture	Up to 0.18	Up to 0.04
Gold	Marker band	Up to 0.14	N/A
Tin-silver solder	Stent	Up to 0.02	N/A
Nitinol	Fenestration reinforcing ring	Up to 0.01	N/A
	Scallop reinforcing	Up to 0.01	N/A
Graphite	Pencil marking	Less than 0.01	Less than 0.01

1.5 Patient population

Adult patients (≥18years) with infrarenal abdominal aortic aneurysms who are not suitable for repair with a standard Zenith Endovascular Graft, when the aneurysmal disease extends up to the level of the renal arteries.

1.6 Intended user

Physicians and teams trained in vascular interventional techniques and in the use of this device.

1.7 Contact with body tissue

The primary device will contact various healthy (non-dissected and non-aneurysmal) and diseased arterial tissues from the femoral/iliac artery to the abdominal aorta. In some cases, contact could be made with a conduit used for vascular access.

During the procedure accessory wire guides and catheters will extend into the aortic arch. It may be necessary to access visceral branch vessels through brachial artery access.

1.8 Operating principle

The Zenith Fenestrated AAA Endovascular Graft is placed over a prepositioned wire guide and advanced to the target anatomy. The radiopaque tip of the sheath of the delivery system is used to visualise the deployment of the stent. The proximal component delivery system is advanced until the radiopaque markers indicating fenestrations and/or scallops on the proximal component are at an appropriate level. The distal component delivery system is advanced until the contralateral limb is above the origin of the contralateral iliac with suitable degree of overlap with the proximal component. Orientation Markers on the graft are checked to ensure the graft is in a satisfactory position prior to deployment (this includes cannulation and releasing the diameter reducing ties of the Zenith® Proximal Component). Stent deployment is controlled by means of release mechanisms and by retraction of the sheath while stabilizing the delivery system.

Post deployment, the fully deployed stent graft components impart an outward radial force upon the inner lumen of the graft, establishing patency in the stented region. The proximal graft incorporates the renal vessels to provide the ability to seal around the infrarenal aorta. When used in conjunction, the proximal and distal graft in combination with the iliac leg provide an intravascular conduit that excludes an abdominal aortic aneurysm from blood flow while maintaining distal perfusion.

1.9 Expected lifetime of device

The graft is not intended to be removed from the patient, therefore the expected lifetime of the device is permanent.

The delivery system is a single use device and has a shelf life of 3 years.

2. INTENDED USE

The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is intended to provide an intravascular conduit that excludes an abdominal aortic aneurysm from blood flow while maintaining distal perfusion.

3. INDICATIONS FOR USE

The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is indicated for the treatment of infrarenal abdominal aortic aneurysms (AAAs) in adult patients (≥18 years) with infrarenal abdominal aortic aneurysms who are not suitable for repair with a standard Zenith Endovascular Graft, when the aneurysmal disease extends up to the level of the renal arteries.

The patient must have aneurysmal morphology suitable for endovascular repair, including:

- Adequate femoral/iliac access compatible with the required delivery systems;
- A non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm with:
 - A length of at least 4mm, i.e. at least 4 mm circumferential sealing zone below the lowest target vessel;
- A diameter measured outer wall to outer wall of no greater than 31mm and no less than 19 mm;
- An angle less than 45 degrees relative to the long axis of the aneurysm; and
- An angle less than 45 degrees relative to the axis of the suprarenal aorta.
- An ipsilateral iliac artery distal fixation site greater than 10 mm in length and 9–21 mm in diameter (measured outer wall to outer wall); and
- A contralateral iliac artery distal fixation site greater than 10 mm in length and 7–21 mm in diameter (measured outer wall to outer wall).

4. CLINICAL BENEFITS

The Zenith Fenestrated AAA Endovascular Graft has been developed to increase survivability of aneurysmal disease and prevent aneurysm rupture.

5. CONTRAINDICATIONS

The Zenith Fenestrated AAA Endovascular Graft is contraindicated in:

- Known allergic reactions to one of the materials in the device should always be considered a contraindication.
- Do not use this device on patients who have a systemic infection or a condition that threatens to infect the endovascular graft, as it may result in serious harm.
- Patients who cannot tolerate contrast agents necessary for intraoperative and post-operative follow-up imaging.
- Patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.
- Patients who cannot tolerate radiation e.g., pregnancy.

6. WARNINGS

6.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Fenestrated grafts are made to a customised design to a specification requested by the responsible Physician, and are tailored to a specific patient's anatomy.
- The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System must only be used by physicians and teams trained in vascular interventional techniques and in the use of this device.
- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Preprocedure imaging reconstruction thickness > 3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenosis from CT.
- **The long-term performance and safety of endovascular grafts has not yet been established.** As a result, life-long, regular follow-up must be undertaken in all patients to assess the ongoing performance of the Zenith Fenestrated AAA Endovascular Graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up.

Specific follow-up guidelines are described in **Section 14, Imaging guidelines and post-operative follow-up.**

- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft.
- At a minimum, annual imaging is required, including:
 - 1) abdominal radiographs to examine device integrity (separation between components, stent fracture or barb separation) and
 - 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease.
- If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.
- The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is not recommended in patients unable to undergo, or who will not be compliant with the necessary preoperative and postoperative imaging and implantation studies as described in **Section 14, Imaging guidelines and post-operative follow-up.**
- Intervention or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak.
- An increase in aneurysm size and/or persistent endoleak may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb and/or leaks may be required to undergo secondary interventions or surgical procedures.
- Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the hemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the gray positioner. After the gray positioner has been removed, if blood loss is excessive, consider placing an uninflated molding balloon or an introduction system dilator within the valve, restricting flow.

6.2 Sterile and single use

- This product is single use and shall not be reused.
- Do not reuse, reprocess or resterilize any part of this device.
- Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death.
- Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another.
- Contamination of the device may lead to injury, illness or death of the patient.

6.3 Patient selection, treatment and follow-up

- The Zenith Fenestrated AAA Endovascular Graft has not been evaluated in the following patient populations:
 - Pregnant or nursing females
- Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease, and/or calcification) should be compatible with vascular access techniques and delivery systems of the profile of a 5.5 mm O.D. (14 French) to 8.6 mm O.D. (22 French) vascular introducer sheath. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolization.
- Key anatomic elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation (> 45 degrees for infrarenal neck to axis of AAA or > 45 degrees for suprarenal neck relative to the immediate infrarenal neck); short proximal aortic neck (< 4 mm); greater than 10% increase in diameter over 15 mm of proximal aortic neck length; and circumferential thrombus and/or calcification at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. Irregular calcification and/or plaque may compromise the fixation and sealing of the implantation sites. Necks exhibiting these key anatomic elements may be more conducive to graft migration.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensable inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.

Considerations regarding the use of the Zenith Fenestrated AAA Endovascular Graft include:

- Risk of aneurysm rupture
- Morbidity and mortality associated with conventional surgical repair
- Comorbidities
- Size of aneurysm
- History of renal failure
- Life expectancy
- Anesthetic risk
- Age of patient
- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 14 French (5.5 mm O.D.) to 22 French (8.6 mm O.D.) vascular introducer sheath.
- Non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm:
 - with a length of at least 4 mm,
 - with a diameter measured outer wall to outer wall of no greater than 31 mm and no less than 19 mm,
 - with an angle less than 45 degrees relative to the long axis of the aneurysm, and
 - with an angle less than 45 degrees relative to the axis of the suprarenal aorta.
- Ipsilateral iliac artery distal fixation site greater than 10 mm in length and 9-21 mm in diameter (measured outer wall to outer wall).
- Contralateral iliac artery distal fixation site greater than 10 mm in length and 7-21 mm in diameter (measured outer wall to outer wall).
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft.

6.4 Implant procedure

- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with 4x4 gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Maintain wire guide position during delivery system insertion.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Fenestrated AAA Endovascular Graft.
- Fluoroscopy should be used during introduction and deployment to confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.
- The use of the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure post-operatively. Care should be taken to limit the amount of contrast media used during the procedure.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).
- Inaccurate placement and/or incomplete sealing of the Zenith Fenestrated AAA Endovascular Graft within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the renal or internal iliac arteries. Renal artery patency must be maintained to prevent/reduce the risk of renal failure and subsequent complications.

Inadequate fixation of the Zenith Fenestrated AAA Endovascular Graft may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.

The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs.

Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system.
- Stop and assess the cause of resistance. Vessel or catheter damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.
- Unless medically indicated, do not deploy the Zenith Fenestrated AAA Endovascular Graft in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Do not cover significant renal or mesenteric arteries (**exception is the inferior mesenteric artery only if the superior mesenteric artery and the anastomosis of Riolan are patent and functioning correctly**) with the endoprosthesis. Vessel occlusion may occur. Clinical studies have identified that with use of the Zenith Fenestrated AAA Endovascular Graft there is a risk of renal impairment/loss (which could lead to a need for dialysis). Clinical studies have shown that this risk is reduced when the renal arteries are stented as part of the procedure to implant a Zenith Fenestrated AAA Endovascular Graft. Stenting of the renal arteries is therefore **strongly** recommended with use of the Zenith Fenestrated AAA Endovascular Graft to reduce the risk of renal impairment/loss.
- Take care during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolisation.
- Before deployment of the suprarenal stent, verify that the position of the access wire guide extends just distal to the aortic arch.
- Verify that the predetermined contralateral iliac leg is selected for insertion on the contralateral side of the patient before implantation.
- Care should be taken not to damage the graft or disturb graft positioning after graft placement in the event reinstrumentation of the graft is necessary.

6.5 Molding balloon use

- Prior to molding in the vicinity of any bridging stent(s) confirm that the aortic section of the stent has been flared.
- Confirm complete deflation of balloon prior to repositioning.
- Do not inflate balloon in iliac vessel outside of graft.

7. PRECAUTIONS

The symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class 1B). However, the device contains cobalt as a stainless steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

7.1 Training requirements for user

CAUTION: Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System **MUST** only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System are outlined below:

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair.
- Knowledge of radiographic image interpretation, device selection and sizing.

A multi-disciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy and repair
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolisation
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimise radiation exposure
- Expertise in necessary patient follow-up modalities

8. MRI SAFETY INFORMATION

 A patient with the Zenith Fenestrated AAA Endovascular Graft may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.	
PARAMETER	NOTES
Item name/identification	Zenith Fenestrated AAA Endovascular Graft
Item manufacturer	Cook Medical
Static magnetic field strength [T]	1.5 T and 3.0 T
Static magnetic field (B ₀) orientation	Horizontal, cylindrical bore
Maximum spatial field gradient [T/m and gauss/cm]	7.2 T/m (720 gauss/cm)
RF excitation	Circularly polarized (CP)
RF transmit coil type	Any transmit coil may be used
Operating mode	Normal operating mode
Maximum whole body SAR [W/kg]	2.0 W/kg
Scan duration	2.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks).
MR image artifact	The presence of the ZFEN stent graft and ancillary devices may produce an image artifact of 68 mm.
If information about a specific parameter is not included, there are no conditions associated with that parameter.	

9. POTENTIAL ADVERSE EVENTS

Use of the Zenith Fenestrated AAA Endovascular Graft poses the following potential risks in addition to those associated with use of a standard Zenith Endovascular Graft:

- Organ impairment/loss due to side-branch vessel occlusion (in particular, renal and/or gastrointestinal impairment/loss).

Clinical studies have identified that with use of the Zenith Fenestrated AAA Endovascular Graft there is a risk of renal impairment/loss (which could lead to a need for dialysis). Clinical studies have shown that this risk is reduced when the renal arteries are stented as part of the procedure to implant a Zenith Fenestrated AAA Endovascular Graft.

Stenting of the renal arteries is therefore **strongly** recommended with use of the Zenith Fenestrated AAA Endovascular Graft to reduce the risk of renal impairment/loss.

Adverse events that may occur and/or require intervention include:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g. buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; barb separation and corrosion
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vessel damage
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)

10. HOW SUPPLIED

STERILE—DO NOT RESTERILIZE—SINGLE USE ONLY.

This device is accompanied by an implant card, that should be given to the patient after it has been completed by the Healthcare Professional. The Zenith Fenestrated AAA Endovascular Graft is supplied sterile (100% ethylene oxide) and pre-loaded in peel-open packages. This product is single use and shall not be reused.

11. INSPECTION OF THE DEVICE

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use the device if damage has occurred, or if both sterile barriers have been damaged or broken. Contact William Cook Australia or your local Cook representative.

If the outer sterile barrier is damaged, inspect the inner barrier prior to use. If the inner barrier is intact, sterile personnel may handle product inside the inner barrier. The inner pouch shall not be handled by sterile personnel or placed on sterile surfaces.

If the inner barrier is damaged but outer barrier is intact, the inner pouch and contents are still considered safe to be handled by sterile personnel and can be placed on sterile surface.

Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient. Do not use after the expiration date printed on the label.

12. DEVICE PREPARATION

Refer to **Section 13.4.1 Proximal body preparation/flush** and **13.4.2 Distal bifurcated body graft preparation/flush**.

13. INSTRUCTIONS FOR USE

Illustrations

13.1 General use information

Prior to use of the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System, review this suggested instructions for use booklet. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary.

These instructions are intended to help guide the physician and do not take the place of physician judgment.

13.2 Pre-implant review

Verify from pre-implant planning that the correct device has been selected. Determinants include:

1. Femoral artery selection for introduction of the main body system, (i.e., define respective contralateral and ipsilateral iliac arteries).
2. Angulation of aortic neck, aneurysm and iliacs.
3. Quality of the aortic neck.
4. Diameters of infrarenal aortic neck and distal iliac vessels.
5. Distance from renal arteries to the aortic bifurcation.
6. Distance from the renal arteries to the hypogastric (internal iliac) arteries/attachment site(s).
7. Aneurysm(s) extending into the iliac arteries may require special consideration in selecting a suitable graft/artery interface site.
8. Consider the degree of vascular calcification.

NOTE: Each respective vessel diameter and length (aorta, ipsilateral iliac and contralateral iliac) provides the necessary criteria for choosing the appropriate endovascular graft.

13.3 Patient preparation

1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
2. Position patient on imaging table allowing fluoroscopic visualisation from the aortic arch to the femoral bifurcations.
3. Expose both common femoral arteries using standard surgical technique.
4. Establish adequate proximal and distal vascular control of both femoral vessels.

13.4 Fenestrated system

13.4.1 Proximal body preparation/flush

1. Remove black-hubbed shipping stylet (from the inner cannula), cannula protector tube (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away® sheath from back of the hemostatic valve. **(Figure 7)** Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove in the tip. **(Figure 8)** Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on extension tube.

NOTE: If there is difficulty in removing the tip protector, brace the device by gripping the gray introducer with thumb and forefinger and apply a gentle side to side twisting motion during removal.

NOTE: Graft flushing solution of Heparinised saline is always used.

NOTE: Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. **(Figure 9)**

2. Attach syringe with normal Heparinised saline to the hub on the inner cannula. Flush until fluid exits the distal tip. **(Figure 10)**

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak 4x4 gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator tip liberally.

13.4.2 Distal bifurcated body graft preparation/flush

1. Remove black-hubbed shipping stylet (from the inner cannula), cannula protector tube (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the hemostatic valve. **(Figure 11)** Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove in the tip. **(Figure 12)** Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on extension tube.

NOTE: If there is difficulty in removing the tip protector, brace the device by gripping the gray introducer with thumb and forefinger and apply a gentle side to side twisting motion during removal

NOTE: Graft flushing solution of Heparinised saline is always used.

NOTE: Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. **(Figure 9)**

2. Attach syringe with normal Heparinised saline to the hub on the inner cannula. Flush until fluid exits the distal tip. **(Figure 13)**

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak 4x4 gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator tip liberally.

13.4.3 Vascular access and angiography

1. Puncture the selected common femoral arteries using standard technique with an 18 or 19UT gage arterial needle. Upon vessel entry, insert:
 - Wire guides – standard .035 inch (0.89 mm) diameter, 145 cm long, J tip or Bentson Wire Guide
 - Appropriate size sheaths (e.g., 6 French (2.0 mm I.D.) or 8 French (2.7 mm I.D.))
 - Flush catheter (often radiopaque sizing catheters – e.g., Centimeter Sizing Catheter or straight flush catheter)
2. Perform angiography to identify level(s) of renals, aortic bifurcation and iliac bifurcations.

NOTE: If fluoroscope angulation is used with an angulated neck it may be necessary to perform angiograms using various projections.

NOTE: A previous planning exercise will have determined which side will be used to introduce the proximal and distal bodies.

13.4.4 Proximal body placement

CAUTION: Verify that the predetermined access site is chosen for the introduction and placement of the proximal body.

1. Ensure the delivery system has been flushed with Heparinised saline and that all air is removed from the system.
2. Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.

NOTE: Monitor the patient's coagulation status throughout the procedure.

3. On ipsilateral side, replace J wire with stiff wire guide (AUS or LES) .035 inch (0.89 mm), 260 cm long and advance through catheter and up to the thoracic aorta. Remove flush catheter and sheath. Maintain wire guide position.

NOTE: A straight angiographic catheter may be inserted up the contralateral side to aid in placement of graft.

4. Before insertion, position proximal body delivery system on patient's abdomen under fluoroscopy to assist with orientation and positioning. Rotate to a position where the anterior markers are situated in the most anterior (12:00 o'clock) position. The sidearm of the hemostatic valve may serve as an external reference to the fenestration(s) and/or scallop(s), anterior and posterior markers and body side markers.

CAUTION: Maintain wire guide position during delivery system insertions.

CAUTION: To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

5. Advance the delivery system until the radiopaque markers indicating the fenestration(s) and/or scallop(s) are at the level of the appropriate arteries. Check that the distal end of the graft is in a satisfactory position above the aortic bifurcation and that the anterior and posterior markers indicate that the graft is in satisfactory orientation. (Figure 14a)
6. The 'tick' marker can be used to assist with orientation of the graft during deployment and prevent 180-degree mis-orientation of the graft.

- A marker formation of ✓ indicates an anterior position of the anterior markers. (Figure 14a, A)
- A marker formation of \ indicates a posterior position of the anterior markers. (Figure 14a, B)

NOTE: Angiography should be performed as needed throughout deployment, to confirm correct placement of the graft.

7. Verify position of the wire guide in the thoracic aorta. Ensure that fenestration(s) and/or scallops are at the level of the appropriate arteries and the anterior markers are in the most anterior (12:00 o'clock) position.

NOTE: The vertical anterior markers, and the horizontal posterior markers should form a cross, on the fluoroscopic image, when correctly oriented. (Figure 14b)

NOTE: The fenestration/scallop markers should be in close apposition to the appropriate side branch vessels.

NOTE: Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. (Figure 9)

Clear identification of fenestration position(s) may not be possible until the graft has been fully unsheathed.

8. Stabilize the gray positioner (the shaft of the delivery system) while withdrawing the sheath. Deploy the first two (2) covered stents by withdrawing the sheath while monitoring device location.
9. Perform angiography and adjust graft placement as necessary. Continue to withdraw the sheath making positional adjustments as necessary.

NOTE: Techniques to ensure that the fenestration(s) and/or scallop(s) will accurately align with their respective vessels will vary, and will depend upon vessel anatomy, graft design, and physician preferences.

10. Proceed with deployment until the graft has been fully unsheathed. (Figure 15)

11. When a satisfactory graft position has been achieved, withdraw the angiographic catheter and wire guide, then exchange to selective wire guide/selective catheter to below the level of the proximal body. Cannulate the partially deployed proximal main body.

NOTE: If a small fenestration is being utilised, care should be taken to properly align the fenestration with the respective vessel.

12. Utilizing contralateral access sheath and wire guide, cannulate and advance a guiding catheter into each small fenestration and its respective vessel. (Figure 16)

NOTE: Non-compliant angioplasty balloons may be used as an alternative to guiding catheters.

NOTE: Cannulation of the scallop and its respective vessel may also be achieved using similar techniques.

NOTE: It is not recommended to use balloons or guiding catheters to guide final placement of large fenestrations as stent struts across fenestration may interfere.

CAUTION: Before release of the diameter reducing ties, verify that the position of the ipsilateral access wire extends just distal to the aortic arch.

CAUTION: During proximal trigger-wire removal, top cap advancement, and subsequent suprarenal stent deployment, verify that the position of the main body wire guide extends just distal to the aortic arch and that support of the system is maximized.

13. Verify proper position of proximal body. Remove the safety lock from the gold trigger-wire release mechanism. Withdraw and remove the trigger-wire to release diameter reducing ties by sliding the gold trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Figure 17)

NOTE: At this point, the proximal main body graft should be fully expanded with the proximal bare stent contained within the top cap.

14. Remove the safety lock from the black trigger-wire release mechanism. Under fluoroscopy, withdraw and remove the trigger-wire to unlock the suprarenal stent from the top cap by sliding the black trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Figure 18)

NOTE: If resistance is felt or system bowing is noticed, the trigger-wire is under tension. Excessive force may cause the graft position to be altered. If excessive resistance or delivery system movement is noted, stop and assess the situation.

If unable to remove the black trigger-wire release mechanism from the top cap, perform the following steps under fluoroscopy:

- a) Remove tension on the trigger-wire by loosening the pin vise and slightly pulling the inner cannula to move the top cap down over the suprarenal stent. Avoid compressing the Zenith Fenestrated proximal body.
- b) Retighten the pin vise.
- c) Remove the black trigger-wire release mechanism.
- d) Continue with (15) in Section 13.4.4 Proximal body placement.

NOTE: If still unable to remove the black trigger-wire release mechanism from the top cap, see Section 13.6 Trigger-wire release troubleshooting.

15. Loosen the pin vise. (Figure 19) Control the position of the graft by stabilizing the gray positioner of the introducer.

CAUTION: Before deployment of suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement and, if required, re-position the access wire guide into the aortic arch to accommodate.

16. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed. (Figures 20a and 20b) Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise to avoid contact with the deployed suprarenal stent.

WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation bars. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

17. Remove the safety lock from the white trigger-wire release mechanism. Withdraw and remove the trigger-wire to detach the distal end of the endovascular graft from the delivery system by sliding the white trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (Figure 21)

NOTE: Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.

13.4.5 Docking of top cap

1. Loosen the pin vise. (Figure 22)
2. Secure sheath and inner cannula to avoid any movement of these components.
3. Advance the gray positioner over the inner cannula until it docks with the top cap. (Figures 23a, 23b and 23c)

NOTE: If resistance occurs, slightly rotate gray positioner and continue to gently advance.

4. Retighten the pin vise and withdraw the entire top cap and gray positioner through the graft and through the sheath by pulling on the inner cannula. (Figure 24) Leave the sheath and wire guide in place.
5. Close the Captor Hemostatic Valve by turning it in a clockwise direction until it stops.

13.4.6 Bridging stent placement and deployment

General use information

In the event that small fenestrations are being utilized, stents may be placed to secure positive alignment.

Standard techniques for placement of bridging stents should be employed during use of stents.

CAUTION: The Zenith Fenestrated endovascular graft is not intended to have a bridging stent placed in its large fenestrations **as stent struts across fenestration may interfere**. The long-term performance and safety of bridging stents in large fenestrations has not been evaluated for this device. Therefore, bridging stents **are not recommended** to be placed within large fenestrations.

1. Return to the guide catheter and wire guide which cannulate the small fenestration and respective vessel.
2. Introduce appropriately sized bridging stent (Section 1.3.4) and advance to the ostium of the fenestration/vessel. Advance into the vessel, leaving approximately 5 mm of stent in the aorta. (Figure 25)

NOTE: Fluoroscopic views tangential to the fenestration will optimize visualization of the stent position relative to the stent graft.

3. Expand stent.

CAUTION: Expansion of the bridging stent should be conducted as directed by the manufacturer.

4. Remove the balloon and replace with an oversized angioplasty balloon. Advance the balloon until the proximal tip is positioned at the ostium.
5. Inflate the balloon to flare the intra-aortic segment of the stent. (Figure 26)

CAUTION: This technique requires high quality imaging. Mobile image intensifiers provide less than adequate imaging quality.

6. Remove the angioplasty balloon.

NOTE: In the event that there is more than one fenestration, repeat the preceding steps for each additional small fenestration.

7. Withdraw renal access sheaths, catheters and wire guides in the contralateral side to a level just above the aortic bifurcation.

13.4.7 Distal bifurcated body placement

1. Ensure the delivery system has been flushed with heparinized saline and that all air is removed from the system.
2. Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.
3. Before insertion, position distal bifurcated body delivery system on patient's abdomen under fluoroscopy to determine the orientation of the contralateral limb. The side arm of the hemostatic valve may serve as an external reference to the contralateral limb radiopaque marker.

NOTE: Distal bifurcated body delivery system will not pass through the sheath used to deliver the proximal body.

NOTE: The proximal body delivery sheath must be removed prior to insertion of the distal bifurcated body delivery system.

4. Insert Distal Bifurcated Body delivery system over the wire, into the femoral artery with attention to sidearm reference.
- CAUTION:** Maintain wire guide position during delivery system insertion.
- CAUTION:** To avoid any twist in the endovascular graft, during any rotation of the delivery, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).
5. Advance delivery system until the contralateral limb is positioned above and anterior to the origin of the contralateral iliac. (**Figure 27**) If the contralateral limb radiopaque marker is not properly aligned, rotate the entire system until it is correctly positioned half way between a lateral and an anterior position on the contralateral side.
6. Repeat angiogram to verify:
 - The degree of overlap with proximal body (no less than 2 stents)
 - The position of the contralateral limb
 - The position of the ipsilateral iliac cuff with respect to the common iliac bifurcation.

Reposition distal bifurcated body as required.

CAUTION: When introducing distal bifurcated body, observe proximal body closely to avoid any disruption to its position.

NOTE: Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. (**Figure 9**)

7. Stabilize the gray positioner (the shaft of the delivery system) while withdrawing the sheath. Deploy the first two (2) covered stents by withdrawing the sheath while monitoring device location. Proceed with deployment until contralateral limb is fully deployed. (**Figure 28**)

NOTE: Tick marker on the contralateral limb of the distal bifurcated body is used to determine anterior/posterior orientation of the contralateral limb. It is not intended to line up with the anterior tick mark on the proximal body.

13.4.8 Contralateral iliac wire guide placement

1. Advance the contralateral catheter and wire guide into the common iliac artery to a level below the short contralateral limb and then manipulate the wire guide into the contralateral limb and into the Distal Bifurcated Body. (**Figure 29**) AP and oblique fluoroscopic views can aid in verification of device cannulation.
2. Advance the angiographic catheter into the body of the graft. Perform angiography to confirm correct position inside the Distal Bifurcated Body. Advance the catheter to where the proximal end of the Distal Bifurcated Body is attached to the introducer.

13.4.9 Distal bifurcated body deployment

1. Perform angiography to confirm proper position of the iliac leg with respect to the internal iliac (hypogastric) artery. Adjust position if necessary.
2. Withdraw sheath until the iliac leg is fully deployed.
3. Remove the safety lock from the black trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the black trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (**Figure 30**) Stop withdrawing sheath.

13.4.10 Iliac leg (contralateral) placement

1. Position the image intensifier to show both the contralateral internal iliac artery and contralateral common iliac artery.
2. Prior to the introduction of the contralateral limb delivery system, inject contrast through the contralateral femoral sheath to locate the contralateral internal iliac artery.
3. Introduce the contralateral iliac leg delivery system into the artery. Advance slowly until the iliac leg graft overlaps at least one full iliac leg stent (i.e., proximal stent of iliac leg graft) inside the contralateral limb of the main body. (**Figure 31**) If there is any tendency for the distal bifurcated body graft to move during this maneuver, hold it in position by stabilizing the gray positioner on the distal bifurcated body component (on the ipsilateral side).

NOTE: If difficulty is encountered advancing the iliac leg delivery system, exchange to a more supportive wire guide. In tortuous vessels the anatomy may alter significantly with the introduction of the rigid wires and sheath systems.

4. Confirm position of distal end of the iliac leg graft. Reposition the iliac leg graft if necessary to ensure both internal iliac patency and a minimum overlap of one full iliac leg stent (i.e., proximal stent of iliac leg graft, maximum overlap of 1.5 stents) within the main body endovascular graft.
5. To deploy, hold the iliac leg graft in position with the gray positioner while withdrawing the sheath. (**Figures 32a and 32b**) Ensure one stent overlap is maintained.
6. Stop withdrawing the sheath as soon as the distal end of the iliac leg graft is released.
7. Under fluoroscopy and after verification of iliac leg graft position, loosen pin vise, retract inner cannula to dock tapered dilator to gray positioner. Tighten pin vise. Maintain sheath position while withdrawing gray positioner with secured inner cannula. (**Figure 33**)
8. Re-check the position of the wire guide.

13.4.11 Distal bifurcated body deployment

1. Remove the safety lock from the white trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the white trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (**Figure 34**)
2. Under fluoroscopy and after verification of iliac leg graft position, loosen pin vise, retract inner cannula to dock tapered dilator to gray positioner. Tighten pin vise. Maintain sheath position while withdrawing gray positioner with secured inner cannula.
3. Re-check the position of the wire guides. Leave sheath and wire guide in place.
4. Close the Captor Hemostatic Valve on the Flexor introducer sheath by turning in a clockwise direction until hemostasis is achieved. (**Figure 35**)

13.4.12 Molding balloon insertion

1. Prepare molding balloon as follows:
 - Flush wire lumen with Heparinised saline.
 - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning it counter-clockwise.
3. Advance the molding balloon over the wire guide and through the Hemostatic Valve of the distal bifurcated body introduction system to level of renal arteries. Maintain proper sheath positioning.

NOTE: Captor Hemostatic Valve may be utilized to assist with hemostasis by turning in a clockwise rotation to the “close” position.

NOTE: Captor Hemostatic Valve should always be in the “open” position when repositioning of molding balloon. (**Figure 9**)

4. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the suprarenal stent and the infrarenal neck, starting proximally and working in the distal direction. (**Figure 36**)

CAUTION: Prior to molding in the vicinity of any bridging stent(s) confirm that the aortic section of the stent has been flared.

CAUTION: Confirm complete deflation of balloon prior to repositioning.

5. Withdraw the molding balloon to the ipsilateral limb distal fixation site and expand.

CAUTION: Do not inflate balloon in iliac vessel outside of graft.

6. Deflate and remove molding balloon. Transfer the molding balloon onto the contralateral wire guide and into the contralateral iliac leg introduction system. Advance molding balloon to the contralateral limb overlap and expand.

CAUTION: Confirm complete deflation of balloon prior to repositioning.

7. Withdraw the molding balloon to the contralateral iliac leg/vessel distal fixation and expand. (**Figure 36**)

CAUTION: Do not inflate balloon in iliac vessel outside of graft.

8. Remove molding balloon and replace it with an angiographic catheter to perform completion angiograms.
9. Remove or replace all stiff wire guides to allow iliac arteries to resume their natural position.

13.5 Final angiogram

1. Position angiographic catheter just above the level of the renal arteries. Perform angiography to verify that the renal arteries are patent and that there are no endoleaks. Verify patency of internal iliac arteries.
2. Confirm there are no endoleaks or kinks and verify position of proximal gold radiopaque markers. Remove the sheaths, wires and catheters.
3. Repair vessels and close in standard surgical fashion

13.6 Trigger-wire release troubleshooting

CAUTION: The following steps should be performed only if unable to remove the proximal trigger-wire as described in **Section 13.4.4 Proximal body placement (14)**.

NOTE: Technical assistance from a Cook product specialist may be obtained by contacting your local Cook representative.

13.6.1 Alternate proximal body deployment

1. Cut the exposed suprarenal stent trigger-wire between the white and black release mechanisms (**Figure 37**) and remove the black release mechanism from the handle.
2. Remove the safety lock from the white (distal) trigger-wire release mechanism.
3. Withdraw the white trigger-wire release mechanism and then remove via its slot over the device inner cannula.

NOTE: This will detach the distal end of the graft from the gray positioner.

4. Using locking forceps, clamp and secure the cut end of the suprarenal stent trigger-wire. (**Figure 38**)
5. Loosen the pin vise and, while maintaining inner cannula and trigger-wire position, advance the gray positioner and sheath into the graft until the tip of the gray positioner is approximately 2 cm from the gold markers on the proximal edge of the proximal body (**Figure 39**). The advanced gray positioner provides added support to the inner cannula.

NOTE: Take care when advancing the gray positioner as there will be sheaths and wire guides positioned in target vessels. Ensure that the tip of the gray positioner is not advanced into the top cap.

NOTE: Maintain gentle tension on the suprarenal stent trigger-wire to remove any slack in the wire as the gray positioner and sheath are being advanced.

6. Lock the pin vise. Confirm that the suprarenal stent trigger-wire is secured by the forceps.
7. Stabilize the gray positioner and slowly advance the sheath until the sheath tip is 2 mm from the gold markers. (**Figure 40**)

NOTE: Take care when advancing the sheath as there will be other sheaths and wire guides positioned in target vessels. Take care not to advance the graft itself during sheath advancement.

8. Stabilize the sheath and slightly retract the gray positioner with inner cannula to move the top cap down over the suprarenal stent. (**Figure 41**)

NOTE: Avoid compressing the body of the graft.

9. Ensure the fenestrations are positioned correctly. Test the resistance on the trigger-wire and make minor adjustments as necessary to reduce trigger-wire resistance (**Section 13.6.1, 8**).
10. Remove the suprarenal stent trigger-wire.
11. Withdraw the sheath until the tapered tip of the gray positioner is exposed.

If there are multiple fenestrations, withdraw all but one of the fenestration guiding catheters. It is the physician's preference as to which guiding catheter to withdraw. It is recommended that the decision be made based on the ease of cannulation of the fenestrations and their respective vessels. (**Figure 42**)

NOTE: Leave the wire guides in place when removing the guiding catheter(s).

12. (a) If there is a single fenestration, then the proximal body needs to be cannulated with a suitable wire guide from the same side as the in-situ guide catheter. This is to allow placement of a molding balloon within the proximal body.

13. Advance a molding balloon along the now available wire guide into the proximal body and position it just superior to the distal-most end of the graft.

NOTE: When using the molding balloon, ensure the appropriate sheath is used – either through the in-situ 20 Fr contralateral sheath or, if direct punctures have been used, through a 14 Fr introducer sheath. This will ensure the safe retrieval of the molding balloon.

14. Inflate the balloon to the full diameter of the graft. (**Figure 43**)

15. Loosen the pin vise (**Figure 44**). Control the position of the graft by stabilizing the gray positioner and balloon catheter.

CAUTION: Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.

16. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed (**Figures 45 and 46**). Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise (**Figure 47**) to avoid contact with the deployed suprarenal stent.

WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

17. If there are multiple fenestrations (**Section 13.6.1, 12**), deflate the balloon then withdraw it leaving the wire guide in place
 - (a) If there is a single fenestration (**Section 13.6.1, 12a**), then the molding balloon and wire guide can be safely removed.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide, which remain in the target vessel

18. Advance the access sheath and aligning stent, which was removed to facilitate the molding balloon, back over the wire guide, through the fenestration and into the respective vessel. (**Figure 48**)

NOTE: Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.

13.6.2 Docking of top cap

1. Loosen the pin vise. (**Figure 49**)
2. Secure sheath and inner cannula to avoid any movement of these components.
3. Advance the gray positioner over the inner cannula until it docks with the top cap. (**Figures 50, 51 and 52**)

NOTE: If resistance occurs, slightly rotate gray positioner and continue to gently advance.

4. Retighten the pin vise and withdraw the entire top cap and gray positioner through the graft and through the sheath by pulling on the inner cannula. (**Figure 53**) Leave the sheath and wire guide in place.
5. Close the Captor Hemostatic Valve on the Flexor introducer sheath by turning it in a clockwise direction until it stops.
6. Return to **Section 13.4.6 Bridging stent placement and deployment**.

13.7 Suprarenal stent deployment troubleshooting

CAUTION: The following steps should be performed only if unable to deploy the suprarenal stent by advancing the top cap as described in **Section 13.4.4, Proximal body placement (15)**.

NOTE: Technical assistance from a Cook product specialist may be obtained by contacting your local Cook representative.

13.7.1 Proximal body placement with distal attachment

If the suprarenal stent cannot be fully deployed by advancing the top cap inner cannula, perform the following steps under fluoroscopy.

1. Tighten the pin-vise. If there are multiple fenestrations, withdraw all but one of the fenestration guiding catheters. It is the physician's preference as to which guiding catheter to withdraw. It is recommended that the decision be made based on the ease of cannulation of the fenestrations and their respective vessels. (**Figure 54**)
 - a) If there is a single fenestration, then the proximal body needs to be cannulated with a suitable wire guide from the same side as the in-situ guide catheter. This is to allow placement of a molding balloon within the proximal body.

NOTE: Leave the wire guides in place when removing the guiding catheter(s).

2. Advance a molding balloon along the now available wire guide into the proximal body and position it just superior to the distal-most end of the graft.

NOTE: When using the molding balloon, ensure the appropriate sheath is used – either through the in-situ 20 Fr Contralateral sheath or, if direct punctures have been used, through a 14 Fr introducer sheath. This will ensure the safe retrieval of the molding balloon.

3. To add support to the inner cannula, inflate the balloon to the full diameter of the graft. (**Figure 55**)
4. Loosen the pin vise. (**Figure 56**)
5. Control the position of the graft by stabilizing the gray positioner and balloon catheter.

CAUTION: Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.

6. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed. (**Figures 57 and 58**) Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise (**Figure 59**) to avoid contact with the deployed suprarenal stent.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s). If the suprarenal stent is fully deployed:

7. a) If there are multiple fenestrations (**Section 13.7.1, 1**), deflate the balloon then withdraw it leaving the wire guide in place. Advance the access sheath and aligning stent, which was removed to facilitate the molding balloon, back over the wire guide, through the fenestration and into the respective vessel. (**Figure 60**)
 - b) If there is a single fenestration (**Section 13.7.1, 1a**), then the molding balloon and wire guide can be safely removed.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

8. Return to **Section 13.6.2 Docking of top cap**.

13.7.2 Proximal body placement without distal attachment

If still unable to fully deploy the suprarenal stent, perform the following steps:

1. Tighten the pin vise and deflate the balloon, while maintaining balloon position.
2. Remove the safety lock from the white trigger-wire release mechanism. Withdraw and remove the trigger-wire to detach the distal end of the endovascular graft from the delivery system by sliding the trigger-wire release mechanism off the handle and remove via its slot over the device inner cannula. (**Figure 61**)
3. Loosen the pin vise (**Figure 62**) and, while maintaining inner cannula position, advance the gray positioner and sheath into the graft until

the tip of the gray positioner is approximately 2 cm from the gold markers on the proximal edge of the proximal body (**Figure 63**). The advanced gray positioner provides added support to the inner cannula.

NOTE: Take care when advancing the gray positioner as there will be sheaths and wire guides positioned in target vessels. Ensure that the tip of the gray positioner is not advanced into the top cap.

4. Lock the pin vise.
5. Verify position of the gold markers and ensure the fenestrations are positioned correctly.
6. To add support to the inner cannula, inflate the balloon to the full diameter of the graft. (**Figure 64**)
7. Loosen the pin vise (**Figure 65**). Control the position of the graft by stabilizing the gray positioner and balloon catheter.

CAUTION: Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.

8. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed (**Figures 66 and 67**). Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise (**Figure 68**) to avoid contact with the deployed suprarenal stent.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

9. a) If there are multiple fenestrations (**Section 13.7.1, 1**), deflate the balloon then withdraw it leaving the wire guide in place. Advance the access sheath and aligning stent, which was removed to facilitate the molding balloon, back over the wire guide, through the fenestration and into the respective vessel. (**Figure 69**)
 - b) If there is a single fenestration (**Section 13.7.1, 1a**), then the molding balloon and wire guide can be safely removed.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

10. Return to **Section 13.6.2 Docking of top cap**.

NOTE: Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.

14. IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

14.1 General

The long-term performance and safety of endovascular grafts has not yet been established. As a result, life-long, regular follow-up must be undertaken in all patients to assess the ongoing performance of the Zenith Fenestrated AAA Endovascular Graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive additional follow-up.

Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs.

Physicians should evaluate patients on an individual basis and prescribe their follow-up relative to the needs and circumstances of each individual patient. The recommended imaging schedule is presented in **Table 14.1**. This schedule continues to be the minimum requirement for patient follow-up and should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.

Annual imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs, non-contrast CT and duplex ultrasound may be used.

- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (separation between components, stent fracture and barb separation).
- Duplex ultrasound imaging may provide information on aneurysm diameter change, endoleak, patency, tortuosity and progressive disease. In this circumstance, a non-contrast CT should be performed to use in conjunction with the ultrasound. Ultrasound may be a less reliable and sensitive diagnostic method compared to CT.

Table 14.1 lists the minimum requirements for imaging follow-up for patients with the Zenith Fenestrated AAA Endovascular Graft. Patients requiring enhanced follow-up should have interim evaluations.

Table 14.1 Recommended imaging schedule for endograft patients

	Angiogram	CT (Contrast and non-contrast)	Abdominal radiographs
Pre-procedure	X ¹	X ¹	
Procedural	X		
Pre-discharge (within 7 days)		X ^{2,3,4}	X
1 month		X ^{2,3,4}	X
3 month		X ^{2,4,5}	
6 month		X ^{2,4}	X
12 month (annually thereafter)		X ^{2,4}	X

¹Imaging should be performed within 6 months before the procedure.

²Duplex ultrasound may be used for those patients experiencing renal failure or who are otherwise unable to undergo contrast enhanced CT scan. With ultrasound, non-contrast CT is still recommended.

³Either pre-discharge or 1 month CT recommended.

⁴If Type I or III endoleak, prompt intervention and additional follow-up post-intervention recommended. See **Section 14.5 Additional surveillance and treatment**.

⁵Recommended if endoleak reported at pre-discharge or 1 month.

14.2 Contrast and non-contrast CT recommendations

- Film sets should include all sequential images at lowest possible slice thickness (≤ 3 mm). DO NOT perform large slice thickness (>3 mm) and/or omit consecutive CT images/film sets, as it prevents precise anatomical and device comparisons over time.
- All images should include a scale for each film/image. Images should be arranged no smaller than 20:1 images on 14 inch x 17 inch sheets if film is used.
- Both non-contrast and contrast runs are required, with matching or corresponding table positions.
- Non-contrast and contrast run slice thickness and interval must match.
- DO NOT change patient orientation or re-landmark patient between non-contrast and contrast runs.

Non-contrast and contrast enhanced baseline and follow-up imaging are important for optimal patient surveillance. It is important to follow acceptable imaging protocols during the CT exam. **Table 14.2** lists examples of acceptable imaging protocols.

Table 14.2 Acceptable imaging protocols

	Non-contrast	Contrast
IV contrast	No	Yes
Acceptable machines	Spiral capable of >40 seconds	Spiral capable of >40 seconds
Injection volume	n/a	150 cc
Injection rate	n/a	>2.5 cc/sec
Injection mode	n/a	Power
Bolus timing	n/a	Test bolus: SmartPrep, C.A.R.E. or equivalent
Coverage - start	Diaphragm	1 cm superior to celiac axis
Coverage - finish	Proximal femur	Profunda femoris origin
Collimation	<3 mm	<3 mm
Reconstruction	2.5 mm throughout – soft algorithm	2.5 mm throughout – soft algorithm
Axial DFOV	32 cm	32 cm
Post-injection runs	None	None

14.3 Abdominal radiographs

The following views are required:

- Four films: supine-frontal (AP), cross-table lateral, 30 degree LPO and 30 degree RPO views centered on umbilicus.
- Record the table-to-film distance and use the same distance at each subsequent examination.

Ensure entire device is captured on each single image format lengthwise. If there is any concern about the device integrity (e.g., kinking, stent breaks, barb separation, relative component migration), it is recommended to use magnified views. The attending physician should evaluate films for device integrity (entire device length including components) using 2-4X magnification visual aid.

14.4 Ultrasound

Ultrasound imaging may be performed in place of contrast CT when patient factors preclude the use of image contrast media. Ultrasound may be paired with non-contrast CT. A complete aortic duplex is to be videotaped for maximum aneurysm diameter, endoleaks, stent patency and stenosis. Included on the videotape should be the following information as outlined below:

- Transverse and longitudinal imaging should be obtained from the level of the proximal aorta demonstrating mesenteric and renal arteries to the iliac bifurcations to determine if endoleaks are present utilizing color flow and color power angiography (if accessible).
- Spectral analysis confirmation should be performed for any suspected endoleaks.
- Transverse and longitudinal imaging of the maximum aneurysm should be obtained.

14.5 Additional surveillance and treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement, 5 mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's comorbidities, life expectancy and the patient's personal choices. Patients should be counseled that subsequent reinterventions including catheter based and open surgical conversion are possible following endograft placement.

15. DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

16. PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

The physician must advise the patient to contact their healthcare provider if he/she experiences changes in their medical condition or in case of concern.

The physician must ensure that the patient has received the implant card, that the patient understands the purpose of the implant card and that the patient is able to access the digital patient facing information. If the patient is not able to access the electronic information, the physician must provide the patient with a printed copy.

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure including:

- Risks and differences between endovascular repair and surgical repair.
- Potential advantages of traditional open surgical repair.
- Potential advantages of endovascular repair.

The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair.

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to post-operative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- The long-term performance and safety of endovascular grafts has not yet been established.** As a result, life-long, regular follow-up must be undertaken in all patients to assess the ongoing performance of the Zenith Fenestrated AAA Endovascular Graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 14 Imaging guidelines and post-operative follow-up**.
- Physicians must advise all patients that it is important to seek prompt medical attention if he/she experiences signs of limb occlusion, aneurysm enlargement or rupture.
- Patients should be counselled on the importance of adhering to the follow-up schedule both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of abdominal aortic aneurysms (AAA). At a minimum, annual imaging and adherence to routine post-operative follow-up requirements is required and should be considered a life-long commitment to the patient's health and well-being.

You can use the basic unique device identifier for this device – 0827002WCA10022021040017X – to find the Summary of Safety and Clinical Performance (SSCP) on the Eudamed website. When the Eudamed is available, use the following link: ec.europa.eu/tools/eudamed

Patient Implant Information

This information will be available for the patient in layman language (except for the MRI information, which is aimed at healthcare professionals) through cookmedical.eu/patient-implant-information. The physician should ensure that the patient has received the implant card giving access to the information. If the patient is not able to access the electronic information, the physician must provide the patient with a printed copy.

17. SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and the competent authority and/or regulatory authority of the country where the device was used.

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© COOK 2026
William A. Cook Australia Pty. Ltd.
Cook Medical Europe Ltd.