

- EN**  
**3**  
**Solus™ Double Pigtail Sof-Flex® Stent and Introduction System**  
 Instructions for Use
- CS**  
**5**  
**Stent Solus™ Sof-Flex® se dvěma pigtaily a zaváděcím systémem**  
 Návod k použití
- DA**  
**7**  
**Dobbeltgrisehale Sof-Flex® stent med Solus™ indføringsystem**  
 Brugsanvisning
- NL**  
**9**  
**Sof-Flex® dubbele pigtailstent met Solus™ introductiesysteem**  
 Gebruiksaanwijzing
- FR**  
**11**  
**Endoprothèse double pigtail Sof-Flex® et système d'introduction Solus™**  
 Mode d'emploi
- DE**  
**13**  
**Solus™-Sof-Flex®-Doppel-Pigtail-Stent mit Einführsystem**  
 Gebrauchsanweisung
- EL**  
**15**  
**Ενδοπρόσθεση Sof-Flex® διπλού σπειροειδούς άκρου Solus™ και σύστημα εισαγωγής**  
 Οδηγίες χρήσης
- IT**  
**17**  
**Stent Solus™ a doppio pigtail in Sof-Flex® con sistema di introduzione**  
 Istruzioni per l'uso
- PT**  
**19**  
**Stent Sof-Flex® de dupla espiral e sistema de introdução Solus™**  
 Instruções de utilização
- ES**  
**21**  
**Stent Sof-Flex® de doble pigtail con sistema de introducción Solus™**  
 Instrucciones de uso
- SV**  
**23**  
**Solus™ Dubbelpigtail Sof-Flex® stent med införingsystem**  
 Bruksanvisning

**INTENDED USE**

This device is used to drain obstructed biliary ducts.

**NOTES**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

**CONTRAINDICATIONS**

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Inability to pass wire guide or stent through obstructed area.

**POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications associated with biliary stent placement include, but are not limited to: trauma to biliary tract or duodenum, obstruction of pancreatic duct, stent migration.

**PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

This stent must be placed under fluoroscopic monitoring.

This stent must only be placed using the Solus Stent Introduction System, which is included with this stent set.

This stent should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Pigtail straightener is used to straighten pigtail curls on stent and to aid introduction of stent into accessory channel. Pigtail straightener is not intended for use in accessory channel.

Assessment must be made to determine necessity of sphincterotomy prior to stent placement. If sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

## WARNING

MRI compatibility of this device has not been established.

## INSTRUCTIONS FOR USE

## Illustrations

1. Advance pigtail straightener, if applicable, along shaft of stent to end. (*See fig. 1*)  
**Note:** Care must be exercised when straightening pigtail curls in order to avoid kinking or breaking the stent.
2. Lubricate guiding catheter with water-soluble lubricant. Back load stent and pigtail straightener onto tip of soft stent introducer. Advance preloaded stent and soft stent introducer over a pre-positioned wire guide.
3. Slide pigtail straightener, if applicable, back along stent, straightening second pigtail curl. Advance stent and remove pigtail straightener as stent enters accessory channel. Dispose of pigtail straightener per institutional guidelines for biohazardous medical waste. **Caution:** Pigtail straightener is not intended for use in accessory channel.
4. Fluoroscopically monitor radiopaque bands on tip of guiding catheter. Radiopaque bands on catheter are arranged in 5 cm increments. When a sufficient length of guiding catheter (a minimum of one radiopaque band) is above majority of stricture, disconnect Luer lock fitting on introducer.  
**Caution:** A sufficient portion of guiding catheter must be above stricture prior to stent positioning.
5. Maintain position of guiding catheter and advance stent into desired position using pushing catheter.
6. Fluoroscopically and endoscopically confirm desired stent placement.
7. Retract and remove both guiding catheter and wire guide into endoscope while maintaining position of pushing catheter.
8. Confirm position of stent once guiding catheter and wire guide are removed.
9. Remove pushing catheter.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

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© COPYRIGHT COOK 2016 ©
**COOK IRELAND LTD.**