
EN 3	Stent Introducer Set Instructions for Use
BG 6	Комплект за въвеждащо устройство за стент Инструкции за употреба
CS 9	Souprava zavadeče stentu Návod k použití
DA 12	Stentindføringsæt Brugsanvisning
DE 15	Stent-Einführset Gebrauchsanweisung
EL 19	Σετ εισαγωγή ενδοπρόσθεσης Οδηγίες χρήσης
ES 22	Equipo introductor de stents Instrucciones de uso
ET 26	Stendi sisestusvahendi komplekt Kasutusjuhend
FR 29	Set d'introducteur d'endoprothèse Mode d'emploi
HR 32	Komplet uvodnice stenta Upute za uporabu
HU 35	Szentbevezető készlet Használati utasítás
IT 39	Set con introduttore per stent Istruzioni per l'uso
LT 42	Stento įstūmimo movos rinkinys Naudojimo instrukcija
LV 45	Stenta ievadītāja komplekts Lietošanas instrukcija
NL 48	Stentintroducerset Gebruiksaanwijzing
NO 51	Stentinnføringsenhetssett Bruksanvisning
PL 54	Zestaw introduktora stentu Instrukcja użycia
PT 58	Conjunto introdutor de stent Instruções de utilização
RO 61	Set de introducere a stentului Instrucțiuni de utilizare
SK 64	Súprava zavádzača stentu Návod na použitie
SL 67	Sklop uvajala stenta Navodila za uporabo
SR 71	Komplet uvodnika stenta Uputstvo za upotrebu
SV 74	Stentinförarsät Bruksanvisning
TR 77	Stent İntrodüser Seti Kullanma Talimatı

EN Figure 1: Stent Introducer Set (SIS)

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TR Şekil 1: Stent Introdüser Seti (SIS)

STENT INTRODUCER SET

CAUTION: U.S Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

This IFU covers Stent Introducer Set (SIS) devices as represented in **Figure 1**. For all stent information, please refer to stent instructions for use. See **Table 1** for compatible Cook biliary stents. The SIS Stent Introducer Set consists of both a guiding catheter and pushing catheter. The guiding catheter is available in 5Fr and 6Fr in a length of 320cm. The pushing catheter is available in 8.5Fr, 10Fr and 11Fr in a length of 170cm.

Performance Characteristics

The function and key features of the device are described below.

Guiding Catheter - The function of the guiding catheter is to guide the stent to its intended location. Guiding catheters contain one or more of the following features:

- Atraumatic tip - The end of the guiding catheter is designed to be atraumatic to the anatomy.
- Radiopacity - The guiding catheter tubing material is radiopaque. The guiding catheter also contains 2 marker bands to further enhance its fluoroscopic visibility.
- Tuohy-Borst - The guiding catheter contains a Tuohy-Borst with a connection to facilitate the attachment of a standard Luer syringe.

Pushing Catheter - The function of the pushing catheter is to advance the stent, over a pre-positioned wire guide and a guiding catheter, to its intended location, and to maintain the position of the stent as it is being deployed. It contains the following feature:

- Proximal Identifier - An identifier at its proximal end. It is used to identify the orientation in which the device must be used.

Device Compatibility

The Stent Introducer Set is compatible with the following:

- Endoscope with 3.2mm, 3.7mm and 4.2mm accessory channel (Refer to label for selection of the appropriate accessory channel)
- 0.035" Wire Guide (Refer to label for selection of the appropriate wire guide size)
- Positioning Sleeve (Flap Protector)
- Contrast Media/Saline
- Standard Luer syringe
- Endoscopic Cap/Wire Guide Locking Device
- Cook Biliary Non-Pigtail Stents (see **Table 1**)

Table 1: Compatible Cook Biliary Non-Pigtail Stents for the Stent Introducer Set

Device Family	Cook Biliary Stents			
	Cotton-Huibregtse®	Soehendra Tannenbaum®	Cotton-Leung®	Plastic Biliary Stent*
	CHBS0	TTS0	CLS0	PBS*
SIS-8.5 (8.5Fr Pushing Catheter & 5Fr Guiding Catheter)	8.5Fr	8.5Fr	8.5Fr	8.5Fr
SIS-10 (10Fr Pushing Catheter & 6Fr Guiding Catheter)	10Fr	10Fr	10Fr	10Fr
SIS-11.5 (11Fr Pushing Catheter & 5Fr Guiding Catheter)	11.5Fr	11.5Fr	11.5Fr	11.5Fr

*Note: PBS devices are not available in all markets, including the US.

Convention is as follows:

- Fr refers to stent diameter in French size.

Patient Population

Adult patients requiring placement of a biliary stent. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with biliary duct obstruction caused by common bile duct stones, malignant biliary obstruction and benign or malignant strictures.

Intended User

The product is intended for use by physicians trained and experienced in ERCP techniques. Standard techniques for the endoscopic placement of stents should be employed.

Contact With Body Tissues

The device is tissue contacting in line with intended use.

Operating Principles

The SIS Stent Introducer Set consists of both a guiding catheter and pushing catheter. The guiding catheter tracks over a pre-positioned wire guide as it guides the stent to its intended location. It has radiopaque marker bands which enhance the fluoroscopic visibility of the guiding catheter. The pushing catheter then advances the stent over the guiding catheter to its intended location while also maintaining the position of the stent as it is being deployed.

INTENDED USE

This device is used for endoscopic biliary stent placement.

INDICATIONS FOR USE

For endoscopic biliary stent placement for biliary drainage of obstructed ducts that could be caused by common bile duct stones, malignant biliary obstruction, benign or malignant strictures or other obstructed biliary conditions requiring drainage.

CLINICAL BENEFITS

The biliary stent introducers have an indirect therapeutic impact as they function to facilitate the successful placement of a biliary stent. The introducers are essential for stent placement.

CONTRAINDICATIONS

- Those specific to ERCP
- Inability to pass a wire guide/guiding catheter/stent through the obstructed area.

WARNINGS

- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

PRECAUTIONS

- A complete diagnostic evaluation must be performed, prior to use, to determine the proper stent size.
- Sphincterotomy is not necessary for the device use.
- Refer to the label for the minimum channel size required for this device.
- Refer to the label regarding the selection of the appropriate wire guide.
- Do not use excessive force to advance the stent.
- The positioning sleeve (supplied with the stent) is not intended for use in the accessory channel of the endoscope.
- The device must be used under fluoroscopic monitoring.
- Dislodgement of a placed stent is possible when attempting additional procedures.
- Do not use this device for any purpose other than stated intended use.
- Store the device in a dry location.
- These introducer sets are only compatible with Cook non-pigtail biliary stents as outlined in **Table 1**.
- Use of this device is restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in ERCP techniques.
- Standard techniques for the endoscopic placement of stents should be employed.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • liver abscess • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary stent placement: fever • obstruction of the pancreatic duct • pain/discomfort • stent migration • trauma to biliary duct or duodenum.

HOW SUPPLIED

These devices are supplied Ethylene Oxide (EO) sterilized in a peel open pouch.

DEVICE PREPARATION

1. Prior to introducing any devices ensure working channel of endoscope is lubricated with water or a water-based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
4. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.
5. Attach endoscopic cap or wire guide locking device to endoscope.

INSTRUCTIONS FOR USE

1. Remove the hub from the end of the guiding catheter, then introduce the guiding catheter into the accessory channel over a pre-positioned wire guide.
2. Advance the guiding catheter until it is positioned above the obstruction. Note: A sufficient length of the guiding catheter must be above the obstruction prior to the stent positioning. Note: The radiopaque band on the distal end of guiding catheter can be used to fluoroscopically determine the catheter position.
3. Load the positioning sleeve onto the duodenal flap end of the stent.
4. Introduce the positioning sleeve and the stent as a unit onto the end of the pre-positioned wire guide and over the guiding catheter. Note: Refer to the Stent IFU for the intended placed orientation of the stent in the biliary duct.
5. Advance the front flap of the stent into the accessory channel, keeping the positioning sleeve over the back flap of the stent until the stent is completely in the accessory channel.
6. Advance the pushing catheter over the pre-positioned wire guide and the guiding catheter. As the pushing catheter advances the stent completely into the accessory channel, slide the positioning sleeve back over the pushing catheter until it reaches the end of the catheter.
7. With the elevator open, advance the pushing catheter in small increments until the stent is in the desired position.
8. Fluoroscopically and endoscopically confirm the desired stent position. Note: If desired, replace the hub on the guiding catheter before flushing with contrast media.
9. After confirming the stent position, gently remove the wire guide and the guiding catheter from the endoscope while maintaining the position of the stent with the pushing catheter.
10. Gently remove the pushing catheter from the accessory channel.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2026
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.