
EN 3	Pushing Catheter Instructions for Use
BG 6	Избутващ катетър Инструкции за употреба
CS 9	Tlačný katetr Návod k použití
DA 12	Skubbekateter Brugsanvisning
DE 15	Platzierungskatheter Gebrauchsanweisung
EL 19	Καθετήρας ώθησης Οδηγίες χρήσης
ES 22	Catéter empujador Instrucciones de uso
ET 26	Suruv kateeter Kasutusjuhend
FR 29	Cathéter pousse-prothèse Mode d'emploi
HR 32	Potisni kateter Upute za uporabu
HU 35	Tolókatéter Használati utasítás
IT 38	Catetere spingitore Istruzioni per l'uso
LT 42	Stumiamasis kateteris Naudojimo instrukcija
LV 45	Bīdītājkatetrs Lietošanas instrukcija
NL 48	Pushing-katheter Gebruiksaanwijzing
NO 51	Skyvekateter Bruksanvisning
PL 54	Cewnik popychający Instrukcja użycia
PT 57	Cateter posicionador Instruções de utilização
RO 61	Cateter de împingere Instrucțiuni de utilizare
SK 64	Zatlačací katéter Návod na použitie
SL 67	Potisni kateter Navodila za uporabo
SR 70	Potisni kateter Uputstvo za upotrebu
SV 73	Påförarkateter Bruksanvisning
TR 76	İtici Kateter Kullanma Talimatı

EN Figure 1: Pushing Catheter (PC)

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TR İtici Kateter

PUSHING CATHETER

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

This IFU covers Pushing Catheter (PC-E) introducers represented by **Figure 1**. See **Table 1** for compatible Cook Biliary Stents. **Table 1** denotes stents only. The Pushing Catheter is an introducer which is 320 cm in length and is available in 5 Fr and 7 Fr which are compatible with stents of sizes 5 Fr and 7 Fr.

Performance Characteristics

The function and key features of the device are described below.

Pushing Catheter - The function of the Pushing Catheter is to advance the stent, over a pre-positioned wire guide to its intended location, and to maintain the position of the stent as it is being deployed. It contains the following feature:

- Proximal Identifier - An identifier at its proximal end. It is used to identify the orientation in which the device must be used.

Endoscopic Visibility – The pushing catheter is visible endoscopically.

Device Compatibility

The Pushing Catheter is compatible with the following:

- Endoscopes (3.2 mm or larger). (Refer to label for selection of the appropriate accessory channel minimum size)
- 0.035" Wire Guides (Refer to label for selection of the appropriate wire guide size)
- Pigtail Straighteners/Positioning Sleeve (Flap Protector)
- Endoscopic Cap/Wire Guide Locking Device
- Cook Biliary Stents (**see Table 1**)

Table 1: Compatible Cook Biliary Stents for the Pushing Catheter

Device Family	Cook Biliary Stents					
	Cotton-Huibregtse®	Cotton-Leung®	Cotton-Leung® (Sof-Flex®)	Plastic Biliary Stent*	Zimmon®	Compass
	CHBSO	CLSO	CLSO-SF	PBS*	ZSO	CBBSO
Pushing Catheter (PC-E)	✓	✓	✓	✓	✓	✓

*Note: PBS devices not available in all markets, including the US.

Convention is as follows:

- Checkmark "✓" means compatible with corresponding stent sizes.

Note: Pushing catheters are compatible with the stents indicated in **Table 1** e.g. PC-7E is compatible with a 7 Fr stent, PC-5E is compatible with a 5 Fr stent etc.

Patient Population

Adult patients requiring placement of a biliary stent. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with biliary duct obstruction caused by common bile duct stones, malignant biliary obstruction and benign or malignant strictures.

Intended User

The product is intended for use by physicians trained and experienced in ERCP techniques. Standard techniques for the endoscopic placement of stents should be employed.

Contact with Body Tissue

This device is tissue contacting in line with intended use.

Operating Principles

The function of the Pushing Catheter is to advance the stent, over a pre-positioned wire guide to its intended location, and to maintain the position of the stent as it is being deployed.

INTENDED USE

This device is used for endoscopic biliary stent placement.

INDICATIONS FOR USE

For endoscopic biliary stent placement for biliary drainage of obstructed ducts that could be caused by common bile duct stones, malignant biliary obstruction, benign or malignant strictures or other obstructed biliary conditions requiring drainage.

CLINICAL BENEFITS

The biliary stent introducers have an indirect therapeutic impact as they function to facilitate the successful placement of a biliary stent. The introducers are essential for stent placement.

CONTRAINDICATIONS

Those specific to ERCP.

Inability to pass wire guide / stent through the obstructed area.

WARNINGS

This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.

Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.

Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

PRECAUTIONS

A complete diagnostic evaluation must be performed, prior to use, to determine the proper stent size.

Refer to the label for the minimum channel size required for this device.

Refer to the label regarding the selection of the appropriate wire guide.

Care must be exercised when straightening the pigtail curls, if any, in order to avoid kinking or breaking the stent.

Do not use excessive force to advance the stent.

The positioning sleeve/pigtail straightener (where supplied with the stent) is not intended for use in the accessory channel of the endoscope.

Dislodgement of a placed stent is possible when attempting additional procedures.

Do not use this device for any purpose other than stated intended use.

Store the device in a dry location.

Use of this device restricted to a trained healthcare professional.

The product is intended for use by physicians trained and experienced in ERCP techniques.

Standard techniques for the endoscopic placement of stents should be employed.

The pushing catheter is only compatible with Cook biliary stents as outlined in **Table 1**.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • liver abscess • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary stent placement: fever • obstruction of the pancreatic duct • pain/discomfort • stent migration • trauma to the biliary duct or duodenum.

HOW SUPPLIED

These devices are supplied Ethylene Oxide (EO) sterilized in a peel open pouch.

DEVICE PREPARATION

1. Prior to introducing any devices ensure working channel of endoscope is lubricated with water or a water-based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
4. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.
5. Attach endoscopic cap or wire guide locking device to endoscope.

Illustrations

INSTRUCTIONS FOR USE

I. FOR NON-PIGTAIL STENTS

1. Gently ensure the full extension of all the flaps.
2. Load the positioning sleeve onto the duodenal flap end of the stent.
3. Introduce the stent, ductal end first, and the positioning sleeve onto a pre-positioned wire guide.
4. Advance the pushing catheter over the wire guide to advance the stent into the accessory channel. As the front flap of the stent enters the accessory channel, keep the positioning sleeve over the back flap until the stent is completely in the accessory channel.

Then refer to steps 4-9 below

II. FOR PIGTAIL STENTS:

1. Straighten the ductal pigtail. Use the pigtail straightener where provided.
2. Introduce the stent, ductal end first, and the pigtail straightener onto a pre-positioned wire guide until the straightener reaches second curl.
3. Advance the pushing catheter over the wire guide to advance the pigtail stent into the accessory channel.
4. As the pushing catheter advances the stent completely into the accessory channel, slide the positioning sleeve/pigtail straightener (where supplied) back over the pushing catheter until it reaches the end of catheter, keeping it clear of the accessory channel.
5. Advance the pushing catheter in small increments until the stent is in the desired position.
6. Fluoroscopically and endoscopically confirm the desired stent position.
7. After confirming the stent position, gently remove the wire guide from the endoscope while maintaining the position of the stent with the pushing catheter.
8. Gently remove the pushing catheter from the accessory channel.
9. These stents may be removed using standard endoscopic techniques.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2026
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.