COOK

MEDICAL

CE
0123

en
3

Dilator

Instructions for Use

bg
6

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Bruksanvisning

DILATOR

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Dilators are designed with a gradual or long distal taper, may include hydrophilic coating, and are available in multiple inner and outer diameter and length combinations. Dilators sizes 13 Fr and larger are made of vinyl radiopaque dilator tubing and dilators smaller than 13 Fr are made of radiopaque polyethylene tubing. Refer to the product label for product specifications (e.g., dilator outer diameter, wire guide endhole size, dilator length).

Performance Characteristics

- Radiopaque dilator material for positioning under imaging techniques
- Tapered dilator tip to facilitate insertion into patient vasculature

Device Compatibility

- Sterile Water
- Heparinized/non-heparinized saline
- Wire guide
- Standard Luer syringe

Patient Population

The target patient population is adult and pediatric patients with a medical condition for which diagnosis or treatment requires dilation of a vascular access site, peripheral vessel, or percutaneous tract as directed by a qualified healthcare provider. The following dilators are designed for the pediatric anatomy and intended for this patient population:

- Hydrophilic coated 12 Fr and smaller
- All non-hydrophilic coated

Intended Users

These devices are used by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices are categorized as external communicating devices. They are in limited contact duration (< 24 hours) with skin, tissue, or circulating blood.

Operating Principle

After initial percutaneous access, dilators are designed to be introduced over a wire guide. The tapered design allows for initial entry; advancing the dilator expands the diameter of the opening and creates a tract. Other devices can then be advanced through the tract. Dilators are designed with tapered open-end distal tip

for a transitional fit with an appropriately sized wire guide and for facilitation of dilatation. Additionally, hydrophilic coating may be applied to the distal portion of the dilator. Some dilators are also offered with a longer tapered tip to make dilating from a smaller french size to a larger french size possible with one dilator. Most dilators have the dilator french size stamped on the proximal fitting; however, dilators manufactured in sizes 20 Fr and larger are designed with shrink tubing attached to the proximal end of the dilator rather than a hub. The french size of the dilator is stamped on the shrink tube.

INTENDED USE

These tapered devices percutaneously dilate vascular puncture sites, peripheral vessels, or catheter tracts in peripheral vascular and non-vascular applications (e.g., biliary, nephrostomy, or gastrostomy tract dilation).

DEVICE AND CLINICAL BENEFITS

These devices dilate peripheral vessel access sites, peripheral vessels, or percutaneous tracts to facilitate completion of subsequent procedures.

INDICATIONS FOR USE

These devices are indicated for adult and pediatric patients with conditions requiring dilation of the peripheral vasculature, skin, and/or subcutaneous tissue to facilitate completion of subsequent percutaneous procedures.

CONTRAINDICATIONS

No contraindications

WARNINGS

- Do not alter the device in any way.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

PRECAUTIONS

- The maximum diameter of the instrument or catheter to be introduced should be considered when selecting the dilator size.
- If resistance is encountered during dilator manipulation, stop and determine the cause before proceeding further.

POTENTIAL ADVERSE EVENTS

- Air embolism (in vascular procedures)
- Bleeding/hematoma
- Infection
- Injury to (or inflammation of) surrounding structures
- Pain

- Subcutaneous emphysema (in non-vascular procedures)
- Thrombus/embolus formation (in vascular procedures)
- Tract leakage (in non-vascular procedures)
- Vascular occlusion (in vascular procedures)
- Vasospasm (in vascular procedures)

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Flush the dilator with sterile water, heparinized or non-heparinized saline (flushing) solution.
2. If the device has hydrophilic coating, immerse dilator in flushing solution to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

NOTE: The surface of the dilator may become dry if not used immediately after activation. Additional wetting with flushing solution will reactivate the hydrophilic coating.

3. Using Seldinger technique, smoothly advance the dilator over a wire guide (not included) to the desired placement within the vessel or tract.
4. Remove dilator from vessel or tract while maintaining the position of the wire guide within the target anatomy for advancement of the next device required for the planned procedure.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2025
Cook Medical Europe Ltd.
William A. Cook Australia Pty. Ltd.