COOK

MEDICAL

CE
0088

EN
3

Wills-Oglesby Percutaneous Gastrostomy Set

Instructions for Use

DA
4

Wills-Oglesby perkutant gastrostomisæt

Brugsanvisning

DE
5

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T _ W P G _ R E V 8

WILLS-OGLESBY PERCUTANEOUS GASTROSTOMY SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Wills-Oglesby Percutaneous Gastrostomy Set consists of a radiopaque polyurethane pigtail catheter with a stiffening cannula for introduction. The catheter has a suture drawstring for catheter loop fixation. The set also contains an appropriately sized catheter needle set, a wire guide, dilators, and accessories for use in percutaneous placement of a catheter for gastroenteric feeding.

INTENDED USE

The Wills-Oglesby Percutaneous Gastrostomy Set is intended for the percutaneous placement of a catheter for gastroenteric feeding.

CONTRAINDICATIONS

None known

WARNINGS

This product contains **NATURAL RUBBER LATEX**, which may cause allergic reactions.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in percutaneous gastroenteric access. Standard techniques for placement of wire guides and catheters should be employed.
- A TFE-coated wire guide must be used with this catheter.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Using a nasogastric catheter, inflate the stomach with air.
2. Using fluoroscopy, choose a catheter insertion site in the left subcostal region that allows the most direct vertical needle path possible. Ideally, the insertion site will be over the rectus sheath and directly over the body of the stomach, toward the greater curvature. **NOTE:** If interposed colon or small bowel anterior to the stomach is suspected, obtain a cross table lateral view.
3. Infiltrate local anesthetic to the level of the peritoneal surface. After intravenous administration of glucagon, maximally distend the stomach with air. **NOTE:** Percutaneous gastropexy may be performed at this point, utilizing suture anchors.
4. Under fluoroscopy, insert the catheter needle set through the abdominal wall and into the stomach, maintaining the most vertical needle position possible.
5. Withdraw the stylet from the needle catheter. **NOTE:** An audible rush of air may occur when the stylet is removed. It is important to avoid decompression of the stomach at this time.
6. Confirm satisfactory needle catheter position by aspirating air or by injecting a small amount of water-soluble contrast medium into the stomach.
7. Insert the wire guide through the needle catheter into the stomach. Typically the wire guide will move along the lesser or greater curvature of the stomach, into the fundus. **NOTE:** The distal tip of the wire guide must be free and mobile rather than lodged against the gastric wall.
8. Remove the needle catheter, leaving the wire guide in place, and use the dilators, in sequence, to dilate the tract.
9. After removing the final dilator, introduce the gastrostomy catheter/stiffening cannula assembly over the wire guide. **NOTE:** This assembly should be partially advanced into the gastric lumen before the catheter can be advanced over the wire guide by itself.
10. After advancing the catheter over the wire guide into the stomach, remove the stiffening cannula and then the wire guide. Confirm intragastric catheter position by injecting a small amount of water-soluble contrast medium, which should outline rugal folds.
11. Form the catheter loop by applying light traction to the loop string. Gently manipulate the catheter until the loop is clearly free within the gastric lumen. **NOTE:** Do not withdraw the loop so far that it begins to enter the gastrostomy tract; doing so may result in withdrawal of catheter sideholes into the tract, and is not necessary for apposition of the stomach to the abdominal wall. **NOTE:** Use fluoroscopy to facilitate proper positioning of the catheter loop.
12. Pull suture end tight to form loop configuration in catheter, and tie securely.
13. Trim off excess suture, and slide latex sleeve over suture to prevent leakage.

14. Attach the included fixation device to the catheter at the skin surface, tape in place, and dress appropriately.
15. Attach an appropriate cap or adapter to the proximal hub of the catheter.
NOTE: Feedings may be initiated the day following catheter placement, after clinical and radiographic evaluations confirm satisfactory catheter position and the absence of complications.
NOTE: Flush catheter with water after each use. If using a Luer tip syringe to flush, remove feeding adapter and connect Luer converter to the hub of the jejunal catheter. Remove Luer converter and connect the feeding adapter after flushing the catheter.

Unlocking Catheter Loop

1. Advance wire guide into catheter.
2. Uncover suture by folding back the latex sleeve, and release suture.
3. Withdraw catheter over wire guide.

NOTE: For catheter exchange, insert new catheter over wire guide. For removal, withdraw wire guide following removal of catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.