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T _ X L C F _ M _ R E V O

EXTRA LARGE CHECK-FLO® INTRODUCER

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Extra Large Check-Flo® Introducers are designed to perform as a guiding sheath and/or introducer sheath. These devices incorporate a radiopaque shaft with a hemostasis valve and are provided with a dilator. These devices are available in various sizes and lengths. Refer to the product label for product specifications (e.g., sheath inner diameter, dilator inner diameter).

Performance Characteristics

- Radiopaque sheath and dilator material for positioning under imaging techniques.
- Hemostasis valve to minimize leakage.
- Rounded sheath tip and tapered dilator tip to facilitate insertion into patient vasculature.

Device Compatibility

- Standard luer syringe
- Heparinized saline
- Contrast media
- 0.035" (0.89 mm) Wire Guide
- Vascular devices with a smaller outer diameter than the specified inner diameter of the sheath

Patient Population

The intended patient population are adult patients with a medical condition whose treatment requires access to the vasculature as directed by a qualified healthcare provider.

Intended User

These devices are intended to be used by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices are for single-use, categorized as external communicating devices in contact with blood for limited (< 24 hours) contact.

Operating Principle

These devices include a sheath and one dilator. The sheaths are available in various lengths. The sheath is made of radiopaque material. The dilator tip extends past the sheath to facilitate introduction and positioning of the introducer assembly. The proximal end of the sheath is comprised of a check-flo hemostasis valve which accepts a range of devices while minimizing blood reflux, as well as a stopcock on the side-arm.

INTENDED USE

The intended purpose of these devices is to serve as an introducer or guiding sheath to establish a pathway for other devices to be placed in the central circulatory system. The introducer is not intended to be placed in the coronary or neuro vasculature.

INDICATIONS FOR USE

These devices are indicated when therapeutic or diagnostic devices are to be introduced into the central circulatory system. The introducer is not intended to be placed in the coronary or neuro vasculature.

CLINICAL BENEFITS

These introducers establish an endovascular pathway to place devices into the central circulatory system (without the introducer itself having a direct diagnostic or therapeutic function).

CONTRAINDICATIONS

- Use during coronary or neurovascular procedures
- Use in pregnant women
- Use in pediatric patients
- Presence of pseudoaneurysm, hematoma, or active infection in the area of the access vessel

WARNINGS

- Before withdrawing the sheath through tortuous anatomy, insert the introducer dilator to avoid possible breakage.
- A large force may be required to remove the dilator from the sheath.

- Do not attempt to insert or withdraw the wire guide and/or introducer if resistance is felt.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- If device use requires the administration of contrast media, adverse reactions may occur. Adverse reactions to contrast media may be mild, moderate, or severe, and they may include allergic-type reactions, adverse physiological reactions, and contrast-induced acute kidney injury. Contrast media should be avoided in patients with acute thyroid storm, and adequate washout time should be allowed in patients who may require radioactive iodine imaging or therapy. A history of adverse reactions to contrast media, and certain comorbidities may increase the risk of exposure to contrast media.
- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

PRECAUTIONS

- This product is intended for use by physicians or healthcare practitioners. Standard techniques for placement of vascular access sheaths should be employed.
- Access vessels should be evaluated to ensure they are adequate for introducer insertion. Factors such as vessel size, tortuosity, or vascular disease may result in kinking or inability to advance the introducer sheath or non-delivery of instruments or catheters to be introduced.
- The maximum diameter of the instrument or catheter to be introduced should be determined to ensure that it will pass through the introducer.
- All instruments or catheters used with this product should move freely through the valve and sheath. Damage to the valve/introducer may result when the fit is tight.
- When inserting, manipulating or withdrawing a device through an introducer always maintain introducer position.
- Before removing or inserting devices through the introducer, aspirate through the side-arm of the valve to clear the introducer, then flush with heparinized saline.
- When inflating a balloon at, or close to, the introducer tip, ensure the balloon is not inside the distal tip of the introducer.
- When puncturing, suturing or incising the tissue near the introducer, use caution to avoid damaging the introducer.
- Ultrasound guided placement of introducers may reduce complications.

POTENTIAL ADVERSE EVENTS

- Access site events (including hematoma, infection, local inflammation, local pain, pseudoaneurysm, seroma)
- Bleeding (including hematoma, hemorrhage [including retroperitoneal hemorrhage], ecchymosis)
- Cardiac rhythm disorders (e.g., AV block, arrhythmias, premature contractions, tachyarrhythmia)
- Compartment syndrome
- Contrast reaction (including contrast-induced nephropathy, allergic reaction, anaphylaxis)
- Death
- Embolism
- Extravasation
- Infection
- Injury to adjacent structures
- Limb ischemia
- Nerve or neurological injury
- Pain
- Stroke or cerebrovascular event
- Subintimal passage
- Thrombosis
- Vasospasm
- Vasovagal reaction
- Vessel injury (including arteriovenous fistula, pseudoaneurysm, vessel dissection, vessel laceration, vessel occlusion, vessel perforation, vessel rupture)

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

Sheath Introduction

1. Upon removal from package, ensure the inner diameter (ID) of the introducer is appropriate for the maximum diameter of the instrument or catheter to be introduced.
2. Using the side-arm of the valve, flush the introducer by filling the introducer assembly completely with heparinized saline.
3. Flush the dilator with heparinized solution.
4. Insert the dilator completely into the introducer, ensuring that the dilator hub is in contact with the introducer hub.
5. Using standard Seldinger technique, access the target vessel with the appropriate needle.
6. Insert an appropriate wire guide into the vessel through the needle, then remove the needle, leaving the wire guide in place.
7. Insert the dilator/sheath combination over the wire guide.
8. Remove the wire guide and dilator; aspirate and flush the introducer side-arm.
9. Insert an appropriately sized device as needed.

Sheath Removal

1. Insert a wire guide at least 10 cm past the tip of the sheath.
2. Insert the introducer dilator over the wire guide into the sheath.
3. Withdraw the sheath and dilator as a unit.
4. Remove the wire guide.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use.

You can use the basic unique device identifier for this device – (0827002CINC050-031_FHX) – to find the Summary of Safety and Clinical Performance (SSCP) on the EUDAMED website. When the EUDAMED is available, use the following link: (ec.europa.eu/tools/eudamed).

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2025
Cook Medical Europe Ltd.
William A. Cook Australia Pty. Ltd.