Zilver® Flex™ 35 Biliary Stents

Instructions for Use

EN
5

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Zilver® Flex™ 35

Mode d'emploi

FR
7

T _ Z I B L L _ R E V . 0

Fig. 1


Fig. 2


Fig. 3


Fig. 4



Fig. 5

ZILVER® FLEX™ 35 BILIARY STENTS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

- a. Handle
- b. Hub
- c. Safety Lock
- d. Introducer Catheter
- e. Introducer Tip
- f. Side-arm Flushing Port
- g. Metal Cannula
- h. Radiopaque Marker on Delivery System
- i. Gold Radiopaque Markers on the Stent

The Zilver Flex 35 Biliary Stents are self-expandable stents made of nitinol. They are flexible, slotted tubes that provide strength and flexibility in the biliary duct upon deployment. Post deployment, the stents impart an outwardly strong radial force upon the inner lumen of the duct, establishing patency in the designated stent region.

The Zilver Flex 35 Biliary Stent is preloaded in a 6.0 French delivery catheter. Stent deployment is controlled by retracting the handle while holding the metal cannula stationary.

MR CONDITIONAL

Non-clinical testing has demonstrated that the Zilver Flex 35 Biliary Stent is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient field of 720 Gauss/cm or less
- Whole body averaged specific absorption rate (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes of scanning (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla), respectively.

In non-clinical testing, the Zilver Biliary Stent produced maximum temperature rises of 0.1, 3.8, 0.8, and 0.1 degrees C (for a single stent at 1.5 Tesla, a pair of overlapping stents at 1.5 Tesla, a single stent at 3 Tesla, and a pair of overlapping stents at 3 Tesla, respectively) at whole body averaged specific absorption rates (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) of MR scanning in a 1.5 Tesla/64 MHz General Electric MR scanner and a 3 Tesla Excite General Electric MR scanner.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Zilver Biliary Stent. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

Heating in the MRI environment for stents with fractured struts is not known.

INTENDED USE

The Zilver Flex 35 Biliary Stents are intended for palliation of malignant neoplasms in the biliary tree.

CONTRAINDICATIONS

- Total biliary obstructions and postoperative strictures that cannot be dilated to permit passage of the introducer catheter
- Stenting of a perforated duct where leakage from the duct could be exacerbated by placement of a stent
- Patients with bleeding disorders
- Severe ascites

WARNINGS

- Possible allergic reactions to nitinol should be considered.
- The safety and effectiveness of this device for use in the vascular system have not been established.

PRECAUTIONS

- These products are intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of percutaneous transhepatic biliary drainage catheters and metal stents should be employed.
- Manipulation of the product requires high-resolution fluoroscopic control.
- The stent has not been designed to inhibit tumor ingrowth. Tissue may grow through stents.
- Do not attempt to remove the stent from the introducer system before use.

- Ensure that the red safety lock is not inadvertently removed until final stent release.
- Deploy the stent over an extra stiff or ultra stiff wire guide.
- Do not attempt to push the handle away from the hub during deployment.
- Do not expose the delivery system to organic solvents (e.g, alcohol).
- Do not use power injection systems with the delivery system.
- Do not rotate any part of the system during deployment.
- The device is intended for single use only. Do not resterilize and/or reuse this device.
- Repositioning of the device after deployment is not possible since the introducer catheter cannot be re-advanced over the stent once deployment begins.
- Placement of this biliary stent requires advanced skills in interventional biliary procedures. The following instructions give technical guidance, but do not obviate formal training in the use of the device.

POTENTIAL ADVERSE EVENTS

Allergic reaction to nitinol

PRODUCT RECOMMENDATIONS

Wire Guide Selection

The use of an extra support .035” wire guide is recommended.

Selection of Stent

Measure the length of the target stricture to determine the length of the stent required. Allow for the area proximal and distal to the tumor to be covered with the stent to protect against interference from further tumor ingrowth.

Measure the diameter of the reference lumen (proximal and distal to the stricture) and refer to Stent Size Selection Table below.

Stent Size Selection Table

Reference Lumen Diameter	Unconstrained Stent Diameter
3.0 – 4.0 mm	5.0 mm
4.0 – 5.0 mm	6.0 mm
5.0 – 6.0 mm	7.0 mm
6.0 – 7.0 mm	8.0 mm
7.0 – 8.0 mm	9.0 mm
8.0 – 9.0 mm	10.0 mm

INSTRUCTIONS FOR USE

Illustrations

1. Determine the proper stent size after complete diagnostic evaluation. The stent deployment must be performed under fluoroscopic control.
2. Introduce the .035” extra support wire guide through the access catheter across the distal segment of the target lesion.
3. Remove the access catheter, leaving the wire guide in place.
4. Use the 1 ml syringe included in the inner package to flush the introducer catheter with saline through the side-arm flushing port.
5. Use the 1 ml syringe to flush the wire guide lumen with saline through the hub.
6. Using fluoroscopic visualization, insert the introducer catheter over the wire guide. Position the catheter so that the distal and proximal radiopaque markers on the stent (i) are at the desired position. The stent is now ready to be deployed. (Fig. 1)

DEPLOYMENT OF THE STENT

1. Before deployment it is important to straighten the proximal part of the introducer catheter as much as possible and to keep the handle in a stable position.
 2. The stent expansion must be performed under fluoroscopic control.
 3. Hold the hub (b) on the metal cannula (g) steady. Prior to deploying the stent, remove the red safety lock (c). (Fig. 2)
 4. Hold the hub end stationary. The stent will be deployed as you pull the handle (a) toward the hub (b). (Figs. 3 and 4)
- CAUTION: The introducer catheter cannot be re-advanced over the stent during deployment.**
5. The stent is fully deployed when the handle (a) reaches the hub (b). (Fig. 5)
 6. Perform a final cholangiogram to ensure proper placement of the device.

NOTE: Stents bridging the papilla should extend no more than 5 mm beyond the papilla and into the duodenum after deployment.
 7. Remove the delivery system. It is recommended to re-advance the outer sheath to its pre-deployment position prior to removing the system.

NOTE: If resistance is met during the withdrawal of the delivery system, re-advance the outer sheath to its pre-deployment position. Withdraw the system as one unit.
 8. After stent placement, a temporary external drainage catheter may be placed over the pre-positioned wire guide. This provides access to the stent if an additional treatment is needed.

Multiple Stent Placement

If placements of multiple stents are required in a patient, to cover the length of the lesion, the following recommendations should be considered:

- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximally to the previously placed stent).
- Stents placed in tandem must overlap to allow for complete coverage of the lesion.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATENTS

This product and its use are protected by one or more of the following patents: United States, 4,580,568; 5,928,280; 5,380,304; 5,700,253; 5,968,088; 6,409,752. Other U.S. patents pending. Foreign patents issued and pending.

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Bloomington, IN 47404 U.S.A.