COOK MEDICAL

CE 0123

EN 3 **Biodesign®**
**Surgisis® Urethral Sling**
Instructions for Use

DE 4 **Biodesign®**
**Surgisis® Urethral-Schlinge**
Gebrauchsanweisung

ES 6 **Biodesign®**
**Cabestrillo uretral Surgisis®**
Instrucciones de uso

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# BIODESIGN® SURGISIS® URETHRAL SLING

## INTENDED USE

The Cook® Biodesign® Surgisis® Urethral Sling is used to provide pubourethral support for the treatment of urinary incontinence resulting from hypermobility or intrinsic sphincter deficiency. The sling is intended for implantation to reinforce soft tissues where weakness exists.

This sling is supplied sterile in peel-open packages and is intended for one-time use.

Rx ONLY This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

SURGISIS® URETHRAL SLING This symbol means the following: Surgisis® Urethral Sling

This product is intended for use by trained medical professionals.

## CONTRAINDICATIONS

This sling is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.

## PRECAUTIONS

- **Do not resterilize.** Discard all open and unused portions of the sling.
- Sling is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard sling if mishandling has caused possible damage or contamination, or if the sling is past its expiration date.
- Ensure that sling is rehydrated prior to securing the sling to the surgical site with sutures, staples, or other fasteners.
- Sling performance has not been evaluated with suture spacing greater than 2 mm.
- Ensure that all layers of the sling are secured when fastening the sling to the surgical site.
- Excess tension on the sling may cause urinary retention in the sling recipient when being treated for urinary incontinence.

## POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials. If any of these conditions occur, the sling should be removed:

- Infection
- Acute or chronic inflammation (initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction

## STORAGE

This sling should be stored in a clean, dry location at room temperature.

## STERILIZATION

This sling has been sterilized with ethylene oxide. Do not resterilize.

## USE OF ANTIMICROBIALS

Because the sling is at times used in surgical fields where sterility cannot be assured, the use of antimicrobials is common practice and may prevent infectious complications.[1] In these cases both antibiotic prophylaxis of the patient and antimicrobial soaking of the sling have been used. Typical flora can be expected to include a variety of aerobic and facultative anaerobic organisms, including, but not limited to, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Pseudomonas aeruginosa*, and *Escherichia coli*. Therefore the following points should be considered:

- Antimicrobials, if used topically or systemically, should provide coverage against a wide spectrum of aerobic and anaerobic organisms.
- Antibacterial prophylaxis, if chosen, should be started prior to surgery and continued post-operatively.[1]

The presence of certain antimicrobials may inhibit revascularization and/or infiltration of cells into the sling. For example, gentamicin is known to hinder neovascularization, epithelialization, and keratinocyte growth,[2-4] while povidone iodine,[5] bacitracin,[3,5] polymyxin B,[6] and vancomycin[7] have all been reported to slow or inhibit wound healing. However, no studies have been conducted to evaluate the combination of antimicrobials with the sling placement.

## INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

**NOTE:** Always handle sling using aseptic technique, minimizing contact with latex gloves.

**REQUIRED MATERIALS**

- A sterile dish (kidney dish or other bowl)
- Sterile forceps
- Rehydration fluid: at least 50 ml of room temperature sterile saline or sterile lactated Ringer's solution for each sling

1. Using aseptic technique, remove the sling inner pouch from its outer bag, and place the inner pouch in the sterile field.
2. Open the inner pouch carefully, and aseptically remove the sling with the sterile forceps.
3. Place the sling into the sterile dish in the sterile field.
4. Add to the dish at least 50 ml of the rehydration fluid for each sling.
5. Allow the sling to rehydrate for at least ten minutes.
6. Prepare the surgical site using standard surgical techniques.
7. Using aseptic technique, trim the sling to fit the surgical site. **NOTE:** An alternative method is to cut the sling to size prior to rehydration. If this method is selected, be sure to rehydrate the sling prior to securing the sling into place. See step 5.
8. Using aseptic technique, transfer the sling to the surgical site and secure the sling into place, avoiding excess tension. **NOTE:** Surgical experience indicates that securing the sling with close tissue approximation produces better outcomes.
9. Complete the surgical procedure.
10. Discard any unused portions of the sling according to institution guidelines for medical waste.

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© COOK BIOTECH 2010
West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.