COOK MEDICAL

CE 0123

en 2 Torcon Angiographic Catheters (Cerebral)
Instructions for Use

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tr 51 Torcon Anjiyografik Kateterler (Serebral)
Kullanma Talimatı

TORCON ANGIOGRAPHIC CATHETERS (CEREBRAL)

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Torcon Angiographic Catheter is a single use, single lumen, braided nylon catheter. These catheters come in a variety of lengths and curves and are compatible with either 0.89 mm (0.035-inch) or 0.97 mm (0.038-inch) wire guides, depending on catheter specification. Refer to the product label for product specifications (e.g., catheter length, distal curve configuration).

Performance Characteristics

- Stainless steel coiling incorporated into shaft material for resistance to kinking.
- Distal tips have an extra soft construction, as well as an atraumatic tapered end.
- Distal curve configurations available to facilitate access into target patient anatomy.
- The Proximal Fitting includes a winged hub and a strain relief.

Device Compatibility

- These catheters are manufactured in 5.0 French and will accept a 0.89 mm (0.035-inch) or a 0.97 mm (0.038-inch) diameter wire guide as well as fit through an appropriately sized sheath.
- Maximum power injection machine limit settings at 1200 psi and product labeled flowrate.
- Compatible with heparinized saline, saline, and contrast media. Contrast media comparable to a viscosity of 11.8 cP at 20°C.

Patient Population

The target population is adult patients who require an angiographic procedure as directed by a qualified healthcare provider.

Intended Users

These devices are intended to be used by physicians or healthcare practitioners (e.g., neurosurgeons, radiologists, or vascular surgeons) with proper authorization to perform cerebral/cervical diagnostic angiography (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices contact the circulatory system for limited duration (≤ 24 hours).

Operating Principle

The catheters are positioned into target vessels via Seldinger technique and serve as a conduit for the injection of contrast agent at the specified flow rate.

INTENDED USE

These catheters facilitate access (via the Seldinger technique) to the cervical and aortic arch arteries and branches or to the veins or venous structures that ultimately connect with the superior vena cava to perform selective catheter-based diagnostic angiography by injection of contrast.

DEVICE AND CLINICAL BENEFITS

The Diagnostic Cerebral Catheters will facilitate angiographic identification of relevant anatomy.

INDICATIONS FOR USE

These devices are indicated for selective catheter-based diagnostic angiography of the cervical and aortic arch arteries/branches or veins/venous structures that ultimately connect with the superior vena cava when an assessment of blood flow or vessel patency of the cerebral/neurovascular space is required.

CONTRAINDICATIONS

No contraindications

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure and/or transmission of disease.
- Do not use the device is the sterile packaging is damaged or unintentionally opened before use.
- Device use requires the administration of contrast media. Adverse reactions to contrast media may be mild, moderate, or severe, and they may include allergic-type reactions, adverse physiological reactions, and contrast-induced acute kidney injury. Contrast media should be avoided

in patients with acute thyroid storm, and adequate washout time should be allowed in patients who may require radioactive iodine imaging or therapy. A history of adverse reactions to contrast media, and certain comorbidities may increase the risk of exposure to contrast media.
- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

PRECAUTIONS

- Manipulation of catheter requires fluoroscopic control.
- Due to thinwall construction, extreme care must be exercised during manipulation and withdrawal. Catheter insertion through a synthetic vascular graft should be avoided whenever possible.
- If resistance is encountered during manipulation, stop and determine the cause before proceeding any further. Prior to advancement of the cerebral catheter, confirm that the vessel can sufficiently accommodate the labeled catheter diameter.
- Do not attempt to heat or reshape the catheter curve; the catheter tip is made from a heat-sensitive material.
- The possible whiplash effect of the long, soft catheter tip must be considered during selective angiography.
- These devices are only for angiographic use; co-axial use and use for endovascular therapeutic interventions are not indicated.

POTENTIAL ADVERSE EVENTS

- Access site-related injury
- Bleeding/hematoma/hemorrhage
- Catheter-related event (e.g., kinking, malfunction)
- Death
- Infection
- Neurological event (e.g., stroke, transient ischemic attack)
- Pain
- Thrombosis/embolism
- Vessel Injury (e.g., perforation, rupture, pseudoaneurysm, dissection, arteriovenous fistula)
- Vessel occlusion
- Vessel spasm

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred. Intended for one-time use.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Inspect for damage and open the inner packaging and device using the standard operating room/catheterization lab aseptic techniques to present the catheter to the operator. Do not use the product if the catheter or packaging is damaged.
2. Flush the device with heparinized saline prior to use.
3. Using standard Seldinger technique, access the target vessel with the appropriate needle.
4. Under fluoroscopic guidance, advance the catheter over an appropriately sized wire guide and/or through an appropriately sized sheath introducer to the intended location. **NOTE:** If resistance is encountered during manipulation, stop and determine the cause before proceeding any further. Prior to advancement of the cerebral catheter, confirm that the vessel can sufficiently accommodate the labeled catheter diameter.
5. When using a contrast power injector, ensure all connections are secure. Please refer to the Product Recommendations section for pressure and flow rate information.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use. You can use the

basic unique device identifier for this device – 0827002CINC012-081_A1Y2 – to find the Summary of Safety and Clinical Performance (SSCP) on the EUDAMED website. When the EUDAMED is available, use the following link: (ec.europa.eu/tools/eudamed).

SERIOUS INCICDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

PRODUCT RECOMMENDATIONS

- Maximum power injector machine settings should not exceed the labeled pressure.
- Labeled flow rate achieved with 11.8cP contrast and contrast temperature of 22° C ± 2° C.
- Select the optimal catheter tip shape and size, considering both the patient's anatomy and the intended site for advancement.

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© COOK 2025
Cook Medical Europe Ltd.
William A. Cook Australia Pty Ltd.