COOK MEDICAL

CE 0123

EN 2 **Flexor®** Introducer
Instructions for Use

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T_FLEXOR-M_REVO

FLEXOR® INTRODUCER

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Flexor® Check-Flo® Introducer Set, Flexor® Check-Flo® Introducer (Ansel Modification), Flexor® Check-Flo® Introducer (Ansel Modification - with High-Flex Dilator and Hydrophilic Coating, Flexor® Check-Flo® Introducer (Ansel 1 Modification – with High-Flex Dilator and Hydrophilic Coating), Flexor® Check-Flo® Introducer (Balkin Up & Over® Contralateral Design), Flexor® Check-Flo® Introducer (Raabe Modification), and Flexor® Tuohy-Borst Side-Arm Introducer (Ansel Modification) devices are designed to perform as a guiding sheath and/or introducer sheath. These devices incorporate a Flexor shaft with a hemostasis valve and are provided with one or more dilators. These devices are available in various sizes, lengths, and configurations. Configurations include differences in shaft construction, such as varying shaft stiffness and distal tip material, dilator material, and hydrophilic coating. Refer to the product label for product specifications (i.e., sheath inner diameter, dilator inner diameter).

Performance Characteristics

- Stainless steel coiling incorporated into sheath material for resistance to kinking.
- Radiopaque marker band incorporated into sheath material for distal tip positioning under imaging techniques.
- Hemostasis valve to minimize leakage.
- Rounded sheath tip and tapered dilator tip to facilitate insertion into patient anatomy.
- If applicable, multiple durometers of nylon sheath shaft material for pushability on proximal end and flexibility on distal end.
- If applicable, hydrophilic coating to enhance lubricity of the device.

Device Compatibility

- Standard Luer syringe
- Heparinized saline
- Contrast media
- 0.46-mm (0.018-inch), 0.89-mm (0.035-inch) and/or 0.97-mm (0.038-inch) diameter wire guide. See product label for specific device compatibility.
- Vascular devices with a smaller outer diameter than the specified inner diameter of the sheath

Patient Population

The target patient population for these devices is patients with a medical condition whose treatment requires access to the vasculature.

Intended Users

These devices are used by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices are for single-use, categorized as external communicating devices in contact with the skin, the vessel wall, the lumen of blood vessels, and the blood, and are for limited (< 24 hours) human contact.

Operating Principle

These devices include a single lumen sheath, and depending on the device specification, one or two dilators. The sheaths are available in various lengths. The sheath is made of a nylon material which is supported by coils to provide flexibility while minimizing kinking. Additionally, hydrophilic coating may be applied to the distal portion of the sheath and dilator. There is a radiopaque marker band located approximately 2-3 mm from the end of the sheath to provide precise positioning. The dilator tip extends past the sheath to facilitate positioning of the introducer assembly. The proximal end of the sheath is comprised of a Check-Flo or Tuohy-Borst haemostasis valve which accepts a range of devices while minimizing blood reflux, as well as a stopcock on the side-arm.

INTENDED USE

The intended purpose of these devices is to serve as a pathway for medical devices to be inserted into the peripheral vascular system.

INDICATIONS FOR USE

These devices are indicated for patients who require medical device placement in the vasculature.

CLINICAL BENEFITS

These devices provide access to the vascular system for the placement of medical devices into the peripheral vasculature.

CONTRAINDICATIONS

- These devices should not be used in the coronary vessels or neurovasculature.
- These devices should not be used in the presence of pseudoaneurysm, hematoma, or infection in the area of the access vessel.

WARNINGS

- If resistance is encountered during advancement of Flexor sheath, assess cause of resistance and consider dilation of any restriction identified or consider alternate treatment strategy. If Flexor sheath is advanced through resistance, force to remove the sheath will be higher, increasing the risk of sheath material or hub separation upon withdrawal.
- Reinsertion of dilator prior to removal of Flexor sheath increases the strength of the sheath and lessens the risk of device separation. If resistance is anticipated or encountered during withdrawal of Flexor sheath, consider carefully reinserting the dilator prior to continuing removal.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- If device use requires the administration of contrast media, adverse reactions may occur. Adverse reactions to contrast media may be mild, moderate, or severe, and they may include allergic-type reactions, adverse physiological reactions, and contrast-induced acute kidney injury. Contrast media should be avoided in patients with acute thyroid storm, and adequate washout time should be allowed in patients who may require radioactive iodine imaging or therapy. A history of adverse reactions to contrast media, and certain comorbidities may increase the risk of exposure to contrast media.

PRECAUTIONS

- This product is intended for use by physicians or healthcare practitioners trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.
- In order to ensure device compatibility, choose a sheath size large enough to accommodate the maximum outer diameter of any devices that will be placed through the sheath.
- All interventional or diagnostic instruments used with this product should move freely through the valve and sheath to avoid damage.
- Do not attempt to heat or reshape the device.
- When inserting, manipulating or withdrawing a device through the sheath, always maintain sheath position.
- When puncturing, suturing or incising the tissue near the sheath, use caution to avoid damaging the sheath.

POTENTIAL ADVERSE EVENTS

- Bleeding
- Compartment syndrome
- Infection
- Injury to surrounding structures
- Pain
- Retroperitoneal bleeding
- Stroke/cerebrovascular event
- Thrombosis/embolism
- Vasospasm
- Vessel injury

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

Sheath Introduction

1. Ensure that the inner diameter (ID) of the sheath is appropriate for the maximum diameter of the instruments to be introduced.
2. Using the side-arm of the valve, flush the sheath by filling the sheath assembly completely with heparinized saline.

3. Flush the dilator with heparinized saline.
4. Insert the dilator completely into the sheath. If the sheath has a Tuohy-Borst valve, tighten the valve around the dilator.
5. If the device has hydrophilic coating, activate the coating by wetting the outer surface of the device with heparinized saline.
 NOTE: For best results, maintain wetted condition of device during placement.
6. Using standard Seldinger technique, access the target vessel with the appropriate needle.
7. Insert a wire guide into the vessel through the needle.
8. Leave the wire guide in place and remove the needle.
9. Insert the dilator/sheath assembly over the wire guide.
10. Remove the dilator, then aspirate and flush through the sheath side-arm.
11. Confirm appropriate sheath position.
12. Insert appropriately sized device(s) as needed.

Sheath Removal

1. Insert a wire guide until its tip extends at least 10 cm past the tip of the sheath.
2. Remove the sheath. Avoid applying traction to hub during removal. If resistance is anticipated or encountered during withdrawal of the Flexor sheath, consider reinserting the dilator and removing the sheath and dilator as a single unit.
3. Remove the wire guide.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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Single Sterile Barrier System with protective packaging inside

EN Single Sterile Barrier System with protective packaging inside • **BG** Единична стерилна бариерна система с вътрешна защитна опаковка • **CS** Systém jednoduché sterilní bariéry s vnitřním ochranným obalem • **DA** Enkelt sterilt barrieresystem med beskyttelsesemballage indvendigt • **DE** Einfaches Sterilbarrieresystem mit innerer Schutzverpackung • **EL** Μονό σύστημα φραγμού αποστείρωσης με προστατευτική εσωτερική συσκευασία • **ES** Sistema de barrera estéril individual con embalaje protector interior • **FI** Yksinkertainen steriilisuojausjärjestelmä, jossa on suojapakkaus sisäpuolella • **FR** Système de barrière stérile simple avec emballage protecteur interne • **HU** Egyszeres steril védőzáras rendszer belső védőcsomagolással • **IT** Sistema a barriera sterile singola con confezionamento protettivo interno • **LT** Viengubo sterilaus barjero sistema su apsaugine vidine pakuote • **LV** Vienas sterilās barjeras sistēma ar aizsargiepakojumu iekšpusē • **NL** Systeem met enkele steriele barrière met daarbinnen een beschermende verpakking • **NO** System med én steril barriere med beskyttende emballasje på innsiden • **PL** System pojedynczej bariery sterylnej z opakowaniem ochronnym wewnątrz • **PT** Sistema de barreira única estéril com embalagem interna protetora • **SK** Systém jednej sterilnej bariéry s vnútorným ochranným obalom • **SL** Enojni sterilni pregradni sistem z notranjo zaščitno embalažo • **SR** Sistem jednostruke sterilne barijere sa zaštitnom unutrašnjom ambalažom • **SV** System med en steril barriär med skyddsförpackning på insidan • **TR** İçinde koruyucu ambalaj olan Tekli Steril Bariyer Sistemi

MD

Medical Device

EN Medical Device • **BG** Медицинско изделие • **CS** Zdravotnický prostředek • **DA** Medicinsk udstyr • **DE** Medizinprodukt • **EL** Ιατροτεχνολογικό προϊόν • **ES** Producto sanitario • **FI** Lääkinnällinen laite • **FR** Dispositif médical • **HU** Orvostechnikai eszköz • **IT** Dispositivo medico • **LT** Medicinos priemonė • **LV** Medicīniskā ierīce • **NL** Medisch hulpmiddel • **NO** Medisinsk utstyr • **PL** Wyrób medyczny • **PT** Dispositivo médico • **SK** Zdravotnícka pomôcka • **SL** Medicinski pripomoček • **SR** Medicinsko sredstvo • **SV** Medicinteknisk produkt • **TR** Tıbbi Cihaz

Australian Sponsor

Australian Sponsor

EN Australian Sponsor • **BG** Възложител за Австралия • **CS** Australský zadavatel • **DA** Australsk sponsor • **DE** Australischer Sponsor • **EL** Χορηγός για την Αυστραλία • **ES** Patrocinador australiano • **FI** Australialainen rahoittaja • **FR** Partenaire australien • **HU** Ausztráliai partner • **IT** Sponsor australiano • **LT** Australijos užsakovas • **LV** Austrālijas sponsors • **NL** Australische opdrachtgever • **NO** Australsk sponsor • **PL** Sponsor australijski • **PT** Promotor australiano • **SK** Zadávateľ z Austrálie • **SL** Avstralski sponzor • **SR** Australijski sponzor • **SV** Australisk sponsor • **TR** Avustralya Sponsoru

EU Importer

EN EU Importer • **BG** Вносител за ЕС • **CS** Dovozce v EU • **DA** EU-importør • **DE** EU-Importeur • **EL** Εισαγωγέας στην ΕΕ • **ES** Importador en la UE • **FI** EU-maahantuoja • **FR** Importateur au sein de l'UE • **HU** EU importőr • **IT** Importatore UE • **LT** Importuotojas į ES • **LV** ES importētājs • **NL** Importeur in de EU • **NO** EU-importør • **PL** Importer UE • **PT** Importador para a UE • **SK** Dovozca do EÚ • **SL** Uvoznik EU • **SR** Uvoznik za EU • **SV** EU-importör • **TR** AB İthalatçısı

Do not use if package is damaged and consult instructions for use

EN Do not use if package is damaged and consult instructions for use • **BG** Да не се използва, ако опаковката е повредена и да се направи консултация с инструкциите за употреба • **CS** Nepoužívejte, pokud je obal poškozený, a prostudujte si návod k použití • **DA** Brug ikke, hvis emballagen er beskadiget. Se brugsanvisningen • **DE** Nicht verwenden, wenn die Verpackung beschädigt ist, und Gebrauchsanweisung beachten • **EL** Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης • **ES** No utilizar si el envase está dañado y consultar las instrucciones de uso • **FI** Ei saa käyttää, jos pakkaus on vahingoittunut; perehdy käyttöohjeisiin • **FR** Ne pas utiliser si l'emballage est endommagé et consulter le mode d'emploi • **HU** Ne használja, ha a csomagolás sérült, és olvassa el a használati utasítást • **IT** Non usare se la confezione è danneggiata e consultare le istruzioni per l'uso • **LT** Nenaudoti, jei pakuotė pažeista, ir žiūrėti naudojimo instrukciją • **LV** Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju • **NL** Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing • **NO** Skal ikke brukes dersom pakningen er skadet; se bruksanvisningen • **PL** Nie używać, jeżeli opakowanie zostało uszkodzone; sprawdzić w instrukcji użycia • **PT** Não utilizar caso a embalagem esteja danificada e consultar as instruções de utilização • **SK** Ak je balenie poškodené, nepoužívajte a prečítajte si návod na použitie • **SL** Če je embalaža poškodovana, pripomočka ne uporabite, glejte navodila za uporabo • **SR** Nemojte koristiti ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu • **SV** Får inte användas om förpackningen är skadad, läs bruksanvisningen • **TR** Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun

cookmedical.com/symbol-glossary

Cook Incorporated
750 Daniels Way
Bloomington, IN 47404
USA

EC REP
Cook Medical Europe Ltd.
O'Halloran Road
National Technology Park
Limerick, Ireland

Australian Sponsor
William A. Cook Australia Pty. Ltd.
95 Brandl St
Eight Mile Plains QLD 4113
Australia

2025-06
T_FLEXOR-M_REV0
cookmedical.com
© COOK 2025