COOK MEDICAL

CE 0123

en 2 **Beacon®** Tip Angiographic Catheters (Visceral)
Instructions for Use

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T_HNBV-M_REV0

BEACON® TIP ANGIOGRAPHIC CATHETERS (VISCERAL)

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Beacon® Tip Angiographic Catheter is a single use, single lumen, braided nylon catheter that contains a radiopaque distal tip. These catheters come in a variety of diameters, lengths, and curves and are compatible with either 0.89 mm (0.035-inch) or 0.97 mm (0.038-inch) wire guides, depending on catheter specification. Some catheters are manufactured with a hydrophilic coating. Refer to the product label for product specifications (e.g., catheter length, distal curve configuration, hydrophilic coating).

Performance Characteristics

- Stainless steel coiling incorporated into shaft material for resistance to kinking.
- Distal tips have a tungsten radiopacifier additive, as well as an atraumatic tapered end.
- Distal curve configurations available to facilitate access into target patient anatomy.
- The proximal fitting includes a winged hub and either a strain relief or cap.
- Some catheters are manufactured with a hydrophilic coating.

Device Compatibility

- These catheters are manufactured in 5.0 French and will accept a 0.89 mm (0.035-inch) or a 0.97 mm (0.038-inch) diameter wire guide as well as fit through an appropriately sized sheath.
- Maximum power injection machine limit settings at 1200 psi and product labeled flowrate.
- Compatible with heparinized saline, saline, and contrast media. Contrast media comparable to a viscosity of 11.8 cP at 20°C.
- 0.89 mm (0.035-inch) labeled catheters: Compatible with 0.89 mm (0.035-inch) embolization coils and appropriately sized microcatheters.
- 0.97 mm (0.038-inch) labeled catheters: Compatible with 0.97 mm (0.038-inch) embolization coils and appropriately sized microcatheters.

Patient Population

The target population is adult patients who require an angiographic procedure as directed by a qualified healthcare provider.

Intended Users

These devices are intended to be used by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices contact the circulatory system for limited duration (≤ 24 hours).

Operating Principle

The catheters are positioned into target vessels via Seldinger technique and serve as a conduit for the injection of contrast agents at the specified flow rate and/or embolization procedures.

INTENDED USE

These catheters facilitate access (via the Seldinger technique) to the visceral and peripheral vessels to perform selective catheter-based diagnostic angiography or venography (by injection of contrast), to introduce microcatheters for super-selective vessel catheterization, or to deploy embolization coils.

INDICATIONS FOR USE

These devices are indicated to facilitate access (via the Seldinger technique) to the visceral and peripheral vessels to perform selective catheter-based diagnostic angiography by injection of contrast, to introduce microcatheters, and/or to deploy embolization coils.

DEVICE AND CLINICAL BENEFITS

These devices facilitate angiographic identification of relevant anatomy, introduce microcatheters, and/or deploy embolization coils.

CONTRAINDICATIONS

None

WARNINGS

- The outer surface of the inner catheter package is **NOT** sterile.

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Device use requires the administration of contrast media. Adverse reactions to contrast media may be mild, moderate, or severe, and they may include allergic-type reactions, adverse physiological reactions, and contrast-induced acute kidney injury. Contrast media should be avoided in patients with acute thyroid storm, and adequate washout time should be allowed in patients who may require radioactive iodine imaging or therapy. A history of adverse reactions to contrast media, and certain comorbidities may increase the risk of exposure to contrast media.
- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

PRECAUTIONS

- These products should be stored in their original sealed foil package with required labeling.
- Manipulation of the catheter requires fluoroscopic control.
- If resistance is encountered during manipulation, stop and determine the cause before proceeding any further. Prior to advancement of the visceral catheter, confirm that the vessel can sufficiently accommodate the labeled catheter diameter.
- Do not attempt to heat or reshape the catheter curve; the catheter tip is made from a heat-sensitive material.
- Activate hydrophilic coating, if present, by wetting the distal end of the catheter with saline. For best results, keep the catheter surface wet during placement.

POTENTIAL ADVERSE EVENTS

- Access site-related injury (e.g., hematoma, vessel occlusion, pseudoaneurysm, arteriovenous fistula)
- Bleeding
- Catheter-related event (e.g., kinking, malfunction)
- Infection
- Pain
- Thrombosis/embolism
- Vessel injury (e.g., dissection, unintended subintimal passage) or injury to surrounding structures
- Vessel spasm

HOW SUPPLIED

This product is sterilized by ethylene oxide gas. Catheters are supplied in a foil pouch and inner peel-open package. The outer surface of the inner package is non-sterile. Keep dry and away from sunlight. Use immediately once foil pouch has been opened. Upon removal from package, inspect the product to ensure no damage has occurred. Intended for one-time use.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Inspect for damage and open foil packaging and remove or expose inner package containing the sterile catheter. Do not use the product if packaging is damaged. **NOTE:** The outer surface of the inner catheter package is not sterile.
2. Inspect for damage and open the inner packaging and device using the standard operating room/catheterization lab aseptic techniques to present the catheter to the operator. Do not use the product if the catheter or packaging is damaged.
3. Activate hydrophilic coating, if present, by wetting the distal end of the catheter with saline. Flush the device with heparinized saline prior to use.
4. Using standard Seldinger technique, access the target vessel with the appropriate needle.
5. Under fluoroscopic guidance, advance the catheter over an appropriately sized wire guide and/or through an appropriately sized sheath introducer to the intended location. **NOTE:** If resistance is encountered during manipulation, stop and determine the cause before proceeding any further. Prior to advancement of the visceral catheter, confirm that the vessel can sufficiently accommodate the labeled catheter diameter.
6. When using a contrast power injector, ensure all connections are secure. Please refer to the Product Recommendations section for pressure and flow rate information.

DISPOSAL OF DEVICES
After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES
These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION
Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use.

SERIOUS INCIDENT REPORTING
If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

PRODUCT RECOMMENDATIONS
- Maximum power injector machine settings should not exceed the labeled pressure.
- Labeled flow rate achieved with 11.9cP contrast and contrast temperature of 22° C ± 2° C.
- Select the optimal catheter tip shape and size, considering both the patient's anatomy and the intended site for advancement.

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© COOK 2025
Cook Medical Europe Ltd.
William A. Cook Australia Pty. Ltd.